

## ROTO SILENTA 630 RS



### **Inhalt des Dokuments / content of the document**

Operating instructions (EN)

Rotoren und Zubehör / Rotors and accessories



# Operating instructions

ROTO SILENTA 630 RS



Translation of the original operating instructions



©2022 – All rights reserved

Andreas Hettich GmbH & Co. KG

Föhrenstrasse 12

D-78532 Tuttlingen, Germany

Telephone: +49 (0)7461 705-0

Fax: +49 (0)7461 705-1125

Email: [info@hettichlab.com](mailto:info@hettichlab.com), [service@hettichlab.com](mailto:service@hettichlab.com)

Internet: [www.hettichlab.com](http://www.hettichlab.com)

## Table of contents

<b>1</b>	<b>About this document.</b>	<b>6</b>
1.1	Use of this document.	6
1.2	Gender reference.	6
1.3	Symbols and labels in this document.	6
<b>2</b>	<b>Safety.</b>	<b>6</b>
2.1	Intended use.	6
2.2	Personnel requirements.	7
2.3	Operator's responsibility.	8
2.4	Safety instructions.	8
<b>3</b>	<b>Device overview.</b>	<b>10</b>
3.1	Technical data.	10
3.2	European registration.	13
3.3	Important labels on the packaging.	14
3.4	Important labels on the device.	14
3.5	Operating and indicator elements.	16
3.5.1	Control.	16
3.5.2	Indicator elements.	16
3.5.3	Controls.	17
3.6	Original spare parts.	19
3.7	Scope of supply.	19
3.8	Returns.	19
<b>4</b>	<b>Transport and storage.</b>	<b>20</b>
4.1	Transport and storage conditions.	20
<b>5</b>	<b>Commissioning.</b>	<b>21</b>
5.1	Unpacking the centrifuge.	21
5.2	Setting up and connecting the centrifuge.	22
5.3	Switching the centrifuge on and off.	23
<b>6</b>	<b>Operation</b>	<b>24</b>
6.1	Opening and closing the lid.	24
6.2	Removing and installing the rotor.	25
6.3	Inserting and removing buckets.	25
6.4	Inserting and removing adapters.	26
6.5	Loading.	26
6.6	Opening and closing the biosafety system.	28
6.6.1	Explanation.	28
6.6.2	Lid with spring-type lock.	29
6.7	Packing instructions, HettLiner.	29
6.8	Centrifugation.	31
6.8.1	Centrifugation in continuous operation.	31
6.8.2	Centrifugation with time preselection.	31
6.8.3	Changing settings during centrifugation.	32
6.9	Quick stop function.	32

<b>7</b>	<b>Software operation</b>	<b>32</b>
7.1	Key switches	32
7.2	Centrifugation parameters	33
7.2.1	Ramp-up and ramp-down parameters	33
7.2.2	Runtime	33
7.2.3	Speed, RPM	34
7.2.4	Integral RCF	34
7.2.5	temperature	34
7.2.6	Relative centrifugal force, RCF	34
7.2.7	Setting the relative centrifugal force (RCF/RZB)	34
7.2.8	Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm <sup>3</sup>	35
7.2.9	Centrifuging radius	35
7.3	Programming	35
7.3.1	Opening or loading programs	35
7.3.2	Entering or changing programs	35
7.3.3	Automatic buffer	36
7.4	Rotor detection	36
7.5	Cooling (for centrifuges with cooling)	36
7.5.1	Instructions, cooling	36
7.5.2	Standby cooling	36
7.5.3	Precooling the rotor	36
7.6	Heating (for centrifuges with heating)	36
7.7	Machine Menu	37
7.7.1	Querying system information	37
7.7.2	Querying operating hours	38
7.7.3	Audible signal	38
7.7.3.1	General	38
7.7.3.2	Setting an audible signal	38
7.7.4	Centrifugation data displayed after switching on	38
7.7.5	Setting the date and time	39
7.8	Program links	39
7.8.1	Linking programs or changing a program link	39
7.8.2	Centrifugation run with program link	40
7.8.3	Deleting program links	41
<b>8</b>	<b>Cleaning and care</b>	<b>41</b>
8.1	Overview table	41
8.2	Cleaning and disinfection instructions	42
8.3	Cleaning	42
8.4	Disinfection	43
8.5	Maintenance	44

<b>9</b>	<b>Troubleshooting</b> .....	<b>45</b>
9.1	Fault description .....	45
9.2	Perform a MAINS RESET .....	46
9.3	Emergency release .....	46
<b>10</b>	<b>Disposal</b> .....	<b>47</b>
10.1	General instructions .....	47
<b>11</b>	<b>Index</b> .....	<b>49</b>

## 1 About this document

### 1.1 Use of this document

- Read this document carefully and in full before commissioning the device for the first time.  
Observe other enclosed instruction sheets where necessary.
- This document is part of the device and must be kept within easy reach.
- This document must be included if the device is passed on to a third party.
- The current version of the document in the available languages can be found on the manufacturer's website: ➔ <https://www.hettichlab.com/de/download-center/>







### 1.2 Gender reference

The employed masculine or feminine language form is to facilitate reading. In the spirit of equal treatment, corresponding terms apply in principle to all genders and do not imply any valuation.

### 1.3 Symbols and labels in this document

#### General symbols

The following markers are used in this document to highlight instructions, results, listings, references and other elements:

Marker	Explanation
1.  2.  3.  ... 	Step-by-step instructions
	Results of action steps
	References to sections of the document and other applicable documents
■ ... ■ ...	Listings without a fixed order
[Buttons]	Controls (for example: buttons, switches)
'Indicator'	Indicator elements (for example: signal lights, screen elements)

## 2 Safety

### 2.1 Intended use

#### Intended use

This device is a laboratory centrifuge suitable for medical applications.

Their exclusive therapeutic purpose is to centrifuge blood in blood bag systems. The separated blood components are transferred by another device (separator) into corresponding satellite bags. The individual components obtained in this way are then used for transfusion or autotransfusion.

The centrifuge is only to be operated by qualified personnel working for blood donation services or hospitals.

The centrifuge is only intended for the uses referred to above.



Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this.

Intended use also includes the observation of all instructions in the user manual and compliance with the required inspection and maintenance intervals.

## Non-intended use

- The centrifuge is not suitable for use in explosive or radioactive, or biologically or chemically-contaminated atmospheres.
- The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.  
The manufacturer generally recommends using only centrifuge tubes with special screw caps for hazardous substances.  
Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- The manufacturer does not recommend centrifugation with flammable or explosive materials.
- The manufacturer does not recommend centrifugation with materials that react chemically with one another with high energy.

## Foreseeable misuse

The manufacturer recommends using only accessories that it has approved for the intended purpose.

Only operate the centrifuge under supervision.

## 2.2 Personnel requirements

### Required qualifications

The user has read the user manual in full and familiarised themselves with the device.



#### NOTICE

##### Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.

### Trained user

The user is trained in laboratory practice and able to carry out the work assigned to them, and to recognise and prevent potential hazards independently.

### Personal protective equipment

Lack of personal protective equipment or unsuitable personal protective equipment increases the risk of impaired health and injury.

- Only use personal protective equipment that is in proper condition.
- Only use personal protective equipment that is adapted to the person (correct size, for example).
- Observe instructions on other protective equipment for specific activities.

## 2.3 Operator's responsibility



*Follow the instructions in this document for proper and safe use of the device.*

*Keep the user manual for future reference.*

### Provide information

- Following the instructions in this document will help:
  - To avoid dangerous situations.
  - To minimise repair costs and downtime.
  - To increase the reliability and service life of the device.
- The operator is responsible for compliance with company regulations, standards and national laws.
- Note and keep the revision of the document separate from the document. If lost, the document can be replaced in the correct revision.
- Keep the user manual available at the place where the device is used.
- Pass the user manual on to the buyer when the device is sold.

### Training of personnel

Lack of knowledge when working with the device may result in serious injury or death.

- Instruct personnel on their tasks and the associated risks in accordance with the instruction.

## 2.4 Safety instructions



### **Reporting serious incidents and notifiable incidents**

*In the event of serious incidents or notifiable incidents involving the device or its accessories, these must be reported to the manufacturer and, where applicable, to the competent authority where the user and/or the patient is registered.*



### **DANGER**

**Risk of contamination for the user due to inadequate cleaning or failure to observe the cleaning instructions.**

- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Observe laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.



### **DANGER**

**Fire and explosion hazard due to hazardous substances in samples.**

- Observe relevant regulations and directives for handling chemicals and hazardous substances.
- Do not use aggressive chemicals (for example: dangerous, corrosive extraction agents such as chloroform, strong acids).



**WARNING**

**Dangers due to insufficient maintenance or maintenance not carried out on time.**

- Follow maintenance intervals.
- Check the device for visible damage or defects.  
If any visible damage or defects are present, take the device out of service and inform a service technician.



**! WARNING**

**Risk of electric shock due to ingress of water or other liquids.**

- Protect the device against external liquids.
- Do not pour any liquids into the interior of the device.
- Transport using original transport packaging.



**! WARNING**

**Contamination with hazardous substances and substance mixtures!**

Observe the following actions for substances and substance mixtures that are toxic, radioactive and/or contaminated with pathogenic microorganisms:

- As a rule, use only centrifuge tubes with special screw caps for hazardous substances.
- Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- If no biosafety system is used, the device is not micro-biologically tight in the sense of standard EN / IEC 61010-2-020.
- Contact the manufacturer if necessary.



**WARNING**

**Risk of injury and damage to the device due to a loose rotor.**

- The driver of the rotor shaft must be correctly seated in the groove of the rotor when mounting the rotor.
- Hand-tighten the nut securing the rotor.
- Check that the rotor is firmly seated.
- Follow maintenance intervals.



**CAUTION**

**Risk of injury due to rotating rotor**

Long hair and items of clothing can get caught on the rotor if the rotor is moved manually.

- Tie long hair back.
- Do not allow garments to hang in the centrifuging chamber.

**NOTICE**

**Damage to the device electronics due to incorrect voltage or frequency at the device circuit breaker.**

- Operate the device with the correct mains voltage and mains frequency.  
The value can be found in the technical data and on the rating plate.

**NOTICE**

**Damage to the device and samples due to premature program termination.**

Premature program termination is caused by power failure, switching off during the program or pulling out the mains plug.

- Do not switch off the device while the program is running.
- Do not trigger the emergency release on the device while the program is running.
- Do not pull out the mains plug while the program is running.

### 3 Device overview

#### 3.1 Technical data

Manufacturer	Andreas Hettich GmbH & Co. KG, D-78532 Tuttlingen		
Model	ROTO SILENTA 630 RS		
Type	5005, 5005-50	5005-80	5005-90
Mains voltage ( $\pm 10\%$ )	400 V 3~ +N		
Mains frequency	50-60 Hz		
power consumption	9700 VA	9400 VA	6600 VA
Power consumption	14 A	13.5 A	9.5 A
Refrigerant	R452A		
max. capacity	12000 ml		
max. permissible density	1.2 kg/dm <sup>3</sup>		
max. speed (RPM)	6000		
max. acceleration (RCF)	6520		
max. kinetic energy	215000 Nm		

Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	yes		
<b>Ambient conditions (EN / IEC 61010-1):</b>			
Installation site	indoors only		
Altitude	up to 2000 m above sea level		
Ambient temperature	5 °C to 40 °C		
Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.		
Overvoltage category (IEC 60364-4-443)	II		
Pollution level	2		
Device protection class	I not suitable for use in potentially explosive atmospheres.		
<b>EMC:</b>			
Emitted EM interference, EM interference immunity	EN / IEC 61326-1 Class B		
Noise level (rotor-dependent)	≤62 dB(A)		≤56 dB(A)
<b>Dimensions:</b>			
Width	813 mm		
Depth	1015 mm	1050 mm	
Altitude	973 mm		
Weight	approx. 355 kg	approx. 367 kg	approx. 306 kg
Manufacturer	Andreas Hettich GmbH & Co. KG, D-78532 Tuttlingen		
Model	ROTO SILENTA 630 RS		
Type	5005-08		
Mains voltage (±10%)	208-220 V +6/-10% 3~ (+N) +PE		
Mains frequency	50-60 Hz		
power consumption	9000 VA		
Power consumption	25 A		

Refrigerant	R452A
max. capacity	12000 ml
max. permissible density	1.2 kg/dm <sup>3</sup>
max. speed (RPM)	6000
max. acceleration (RCF)	6498
max. kinetic energy	215000 Nm
Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	yes
<b>Ambient conditions (EN / IEC 61010-1):</b>	
Installation site	indoors only
Altitude	up to 2000 m above sea level
Ambient temperature	5 °C to 30 °C
Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.
Overvoltage category (IEC 60364-4-443)	II
Pollution level	2
Device protection class	I not suitable for use in potentially explosive atmospheres.
<b>EMC:</b>	
Emitted EM interference, EM interference immunity	EN / IEC 61326-1 Class B
Noise level (rotor-dependent)	≤62 dB(A)
<b>Dimensions:</b>	
Width	813 mm
Depth	1015 mm
Altitude	973 mm
Weight	approx. 401 kg

## Rating plate




Fig. 1: Rating plate

- 1 Item number
- 2 Serial number
- 3 Revision
- 4 Equipment number
- 5 Data matrix code
- 6 any labelling indicating whether medical device or in vitro diagnostic medical device
- 7 Global Trade Item Number (GTIN)
- 8 Date of manufacture
- 9 Serial number
- 10 any EAC mark, CE mark
- 11 Country of manufacture
- 12 Date of manufacture
- 13 Mains frequency
- 14 Maximum kinetic energy
- 15 Maximum permissible density
- 16 Manufacturer's address
- 17 any Coolant circuit pressure
- 18 any Coolant capacity
- 19 any Coolant type
- 20 Revs per minute
- 21 Performance values
- 22 Mains voltage
- 23 any Device designation
- 24 Manufacturer's logo

### 3.2 European registration

#### Device conformity



Device conformity according to EU directives.

#### Notified body:

mdc medical device certification GmbH – Notified Body CE 0483

Tel: +49 (0)711 253597 0

Fax: +49 (0)711 258597 10

E-mail: [mdc@mdc-ce.de](mailto:mdc@mdc-ce.de)

Website: [www.mdc-ce.de](http://www.mdc-ce.de)

Address: Kriegerstrasse 6, D-70191 Stuttgart, Germany


#### Single Registration Number

SRN: DE-MF-000010680


Basic-UDI-DI

Basic-UDI-DI	Device assignment
040506740100019J	ROTO SILENTA 630 RS (medical device)


### 3.3 Important labels on the packaging



**TOP**  
This is the correct upright position of the shipping container for transport and/or storage.




**FRAGILE GOODS**  
The contents of the shipping container are fragile, so it must be handled with care.



**PROTECT FROM MOISTURE**  
The shipping container must be kept away from rain and kept in dry conditions.




**TEMPERATURE LIMITATION**  
The shipping container must be stored, transported and handled within the indicated temperature range (-20 °C to +60 °C).




**HUMIDITY LIMITATION**  
The shipping container must be stored, transported and handled within the indicated air humidity range (10% to 80%).

nicht kondensierend  
non-condensing  
sans condensation



**STACK LIMITATION BASED ON QUANTITY**  
Maximum number of identical packages that may be stacked on the lowest package, "n" standing for the number of packages allowed. The lowest package is not included in "n".

### 3.4 Important labels on the device



*The labels on the device must not be removed or covered, or have anything pasted over them.*





Attention, general danger area.

Ensure you read the instructions for commissioning and operation and observe the safety instructions before using the device.



Biohazard warning.




Warning: hot surface.

Failure to observe this warning may result in damage to property and/or personal injury.




Direction of rotation of the rotor.

The orientation of the arrow indicates the rotor's direction of rotation.



Symbol for the separate collection of electrical and electronic equipment, in accordance with Directive 2012/19/EU (WEEE).

Use in European Union countries, Norway and Switzerland.



Key switch positions.



The centrifuge is equipped with an optical interface.

The optical interface is marked with a symbol.

The centrifuge can be controlled and data retrieved via the interface. The *[PROG]* button lights up during data communication.



Equipotential: Connector (PE plug) for equipotential bonding (only for centrifuges with a PE plug).

### 3.5 Operating and indicator elements

#### 3.5.1 Control




Fig. 2: Control

#### 3.5.2 Indicator elements



Fig. 3: [Lid] button

- The button lights up when the lid is closed.



Fig. 4: 'Lid closed' indicator

- The indicator appears when the lid is closed.



Fig. 5: 'Lid open' indicator

- The indicator appears when the lid is open.

#### LOCK 1, LOCK 2

Fig. 6: [Key switch position] indicator

- The indicator appears when the key switch is in this switch position.

#### LOCK 4, LOCK 5

Fig. 7: [Key switch position] indicator

- The indicator appears when the program lock is enabled during serial communication (for centrifuges with serial communication only).

#### PC, ~~PC~~

Fig. 8: [Serial communication] indicator

- The indicator appears if the centrifuge has a serial interface and the centrifuge is connected respectively not connected.



Fig. 9: 'Rotation' indicator

- The indicator appears when the rotor is turning.

## STOP

Fig. 10: [STOP] indicator

- The indicator appears during the centrifugation run as long as the rotor is turning.  
The indicator flashes after an emergency stop.
- The indicator flashes after an emergency stop.

### 3.5.3 Controls



Fig. 11: [Rotary knob]

- Setting the individual parameters.  
Turning anticlockwise decreases the value.  
Turning clockwise increases the value.



Fig. 12: [Mains switch]

- Switch the device on and off.




Fig. 13: [Key switch]

- The key switch switches various functions on and off, depending on the position.




Fig. 14: [Temperature and centrifuging radius] button

- Temperature setpoint, parameter T/°C  
Adjustable from -20°C to +40°C, in 1°C increments (adjustable from -20°C to +90°C with heating/cooling option).  
The lowest achievable temperature is rotor dependent.

- Centrifuging radius  
Parameter r/mm. Input in mm.



Fig. 15: [Ramp-up parameters] button

- Ramp-up levels, parameters  
Level 9 = shortest ramp-up time, Level 1 = longest ramp-up time.
- Ramp-up time parameters  
The adjustable time range is dependent on the set speed.



Fig. 16: [Ramp-down parameters] button

- Brake levels, parameters  
R = Linear braking curve,  
B = similar to an exponential braking curve.  
Level R9, B9 = short ramp-down time, ...  
Level R1, B1 = long ramp-down time,  
Level R0 = unbraked ramp-down.
- Ramp-down time, parameters  
The adjustable time range is dependent on the set speed.
- Brake cut-off speed, parameter n<sup>(\*)</sup> /RPM  
Unbraked ramp-down takes place after reaching this speed.



Fig. 17: [Lid] button

- Open the lid.



Fig. 18: [RCF] button

- Querying the integral RCF, parameter ∫RCF




Fig. 19: [n] button

- Speed, parameter RPM.  
Adjustable from 50 RPM to the maximum rotor speed (n-max-Rotor)
- Querying the maximum rotor speed, parameter n-max-Rotor



Fig. 20: [PROG] button

- Select program location, parameter PROG No.



Fig. 21: [RCF] button

- Relative centrifugal force, parameter RCF/RZB.  
A numerical value can be set that gives a speed between 50 RPM and the maximum rotor speed (n-max-Rotor). Adjustable in 1 second increments.
- Querying the maximum RCF of the rotor, parameter RCF-max-Rotor.



Fig. 22: [RCL] button

- Retrieving programs.



Fig. 23: [START] button

- Start centrifugation run.
- Acceptance of changes during the centrifugation run.



Fig. 24: [STO] button

- Saving programs. 89 programs can be saved (program locations 1 to 89).  
The program locations "----" and 90 to 99 serve as an automatic buffer.  
No programs can be stored in these program locations.



Fig. 25: [STOP] button



Fig. 26: [t] button

- End the centrifugation run.  
The rotor ramps down to a stop at the preselected brake level.

- Runtime, parameter t/min:sec  
Parameter t/min: Adjustable from 1 to 999 min, in 1 minute increments.  
Parameter t/ :sec Adjustable from 1 - 59 s, in 1 second increments.  
Continuous operation "---:--"

### 3.6 Original spare parts

Only use original spare parts from the manufacturer and approved accessories.

### 3.7 Scope of supply

The following accessories are supplied with the centrifuge:

- 1 grease for the trunnions
- 1 open-ended spanner (17 mm and 19 mm AF)
- 1 Hex key (5 mm x 170)
- 10 covering caps Ø12
  
- 3 Wood screw
- 3 Washer
- 2 Metal rail
- 4 Roofing nail
- 1 user manual
- 1 Instruction sheet for packaging removal
- 1 Instruction sheet for installation
- 3 Program data sheet for S control unit

Additionally for types 5005-08, 5005-80, 5005-90:

- 1 Notes on setup and installation

Additionally for delivery in Germany:

- 1 inspection book

Rotors and the corresponding accessories are supplied depending on the order.

### 3.8 Returns

An original Return Material Authorisation (RMA) form from the manufacturer must always be requested for a return. Secure and reliable acceptance and booking in of the goods with the manufacturer is not possible without an original RMA form from the manufacturer. The Return Material Authorisation (RMA) form contains a Declaration of No Objection (UBE), which must be completed in full and enclosed with the return.

If the device and/or accessories are returned to the manufacturer, the complete return shipment must be cleaned and decontaminated by the sender. If returns are not cleaned and/or decontaminated or are insufficiently cleaned and/or decontaminated, this will be performed by the manufacturer and charged to the sender.

The original transport locks must be attached for return shipment, see → *Chapter 4 'Transport and storage' on page 20*. The device must be shipped in its original packaging.

## 4 Transport and storage

### 4.1 Transport and storage conditions

#### Transport conditions



#### NOTICE

**Damage to the device due to failure to use the transport locks.**

- Secure the transport locks before transporting the device.



#### NOTICE

**Damage to the device due to condensation.**

There is a risk of condensation forming on electrical components when component surfaces are cold and the surrounding air is warmer. The condensation that forms may cause a short circuit and/or destroy electronics.

- Warm the device up for at least 3 hours in a warm room before connecting it to the mains.  
or
- Warm up for 30 minutes in a cold room.

- Before transporting, fasten the transport lock and disconnect the device from the mains socket.
- The transport temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.
- Be aware of the weight of the device.
- When transporting using a transport aid (e.g., a pallet truck), the transport aid must be able to carry at least 1.6 times the transport weight of the device.
- Secure the device to prevent it tipping over and falling down during transport.
- Never transport the device sideways or upside down.

#### Storage conditions

- The device must be stored in the original packaging.
- Only store the device in dry rooms.
- The storage temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.

## 5 Commissioning

### 5.1 Unpacking the centrifuge



#### CAUTION

Danger of crushing due to parts falling out of the transport packaging.

- Keep the device balanced during the unpacking process.
- Only open the packaging at the points provided for this purpose.



#### CAUTION

Risk of injury from lifting heavy loads.

- Provide an adequate number of helpers.
- Note the weight. See → Chapter 3 'Device overview' on page 10.



#### NOTICE

Damage to the device due to improper lifting.

- Do not lift the centrifuge by the control panel or the control panel holder.

#### Personnel:

- Trained user







Fig. 27: Unpacking

- 1 Flat
- 2 Device feet
- 3 Metal rail
- 4 Wooden beams

1. ➤ Remove the packaging.
2. ➤ Remove the wooden beam (4).
3. ➤ Attach the metal rails (3) to the wooden pallet using two nails each.
4. ➤ Slide the wooden beam (4) under the metal rails (3) to support them.
5. ➤ Place an open-end spanner (size 10 mm) on the surfaces (1) and turn the device feet (2) upwards as far as possible.
6. ➤ Carefully roll the centrifuge off the wooden pallet over the metal rails (3).

7.  Push the centrifuge to its installation location.
8.  Place the open-end spanner (size 10 mm) on the surfaces (1) and turn the feet (2) down until the castors are no longer in contact with the ground.
9.  Align the centrifuge so it is horizontal by turning the device feet (2).

## 5.2 Setting up and connecting the centrifuge

### Setting up the centrifuge



#### WARNING

Risk of injury due to failing to maintain a sufficient distance to the centrifuge.

- As per EN / IEC 61010-2-020, no persons, hazardous materials or objects may be present within a **safety zone of 300 mm** around the centrifuge during a centrifugation run.
- A distance of **300 mm** from the ventilation slots and ventilation openings of the centrifuge must be maintained.



#### CAUTION

Risk of crushing and damage to the device due to it falling down because of vibration-induced position alterations.

- Place the device on a stable and level surface.
- Select the installation surface dependent on the weight of the device.






#### NOTICE

Damage to the samples and the device if the ambient temperature exceeds or falls below the respective maximum/minimum permissible ambient temperature.

- Comply with the maximum and minimum permissible ambient temperatures for installation of the device.
- Do not place the device next to a heat source.
- Do not expose the device to direct sunlight.
- Do not expose the device to frost.

#### Personnel:

- Trained user

1.  Place the device on a stable and level surface.
2.  Maintain a distance of 300 mm around the device.
3.  Comply with the ambient conditions in the technical data (→ *Chapter 3.1 'Technical data' on page 10*).

### Connecting the centrifuge



#### NOTICE

Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.



**NOTICE****Damage to the device due to condensation.**

There is a risk of condensation forming on electrical components when component surfaces are cold and the surrounding air is warmer. The condensation that forms may cause a short circuit and/or destroy electronics.

- Warm the device up for at least 3 hours in a warm room before connecting it to the mains.  
or
- Warm up for 30 minutes in a cold room.

**Personnel:**

- Trained user

1. ▶ Types 5005-90 and 5005-08 are permanently connected devices.  
With permanently connected devices, a switch must be fitted in the building installation to disconnect the mains supply to the device, in accordance with the laboratory equipment standard EN / IEC 61010-1.  
  
The switch must be located near the device, be readily accessible to the user and be marked as the means of isolation for this device.  
It must be possible to secure the switch to prevent it being switched on again.
2. ▶ A type B residual current circuit breaker must be used if the device is additionally protected with a residual current circuit breaker in the building installation.  
  
When using a different type, the residual current circuit breaker may either not switch off the unit if there is a fault on the unit, or it may switch off the unit even though there is no fault on the unit.
3. ▶ Types 5005-08, 5005-80 and 5005-90 must be connected in accordance with the instructions for setup and installation (AH5005-02).
4. ▶ Centrifuge with PE connector:  
If necessary, connect the PE connector on the back of the device to an additional medical potential equalisation system.
5. ▶ Check whether the mains voltage matches the specification on the rating plate.
6. ▶ For types 5005, 5005-50 and 5005-80:  
Connect the device to a standard mains socket using the mains cable.

### 5.3 Switching the centrifuge on and off.

**Switching the centrifuge on****Personnel:**

- Trained user

- Set the mains switch to *///*.
  - ➔ The buttons flash, depending on the centrifuge type.

The following indicators appear one after the other, depending on the centrifuge type:

  - the centrifuge model
  - the last rotor code recognised by the rotor detection and the maximum rotor speed
  - the program version
  - When the lid is closed: 'OPEN OEFFNEN' indicator
  - When the lid is open: The centrifugation data of the last program used or program 1.

**Immediate display of centrifugation data after switching on**

1. → Set the mains switch to *///*.
2. → Press any button (except the *[STOP]* button) at the first visual change in the display (inverse display).
  - ➔ Centrifugation data is displayed.

**Switching off the centrifuge**

- The rotor is stationary.
- Set the mains switch to *[0]*.

## 6 Operation

### 6.1 Opening and closing the lid

**Opening the lid**

- Personnel:**
- Trained user
- The centrifuge is switched on.  
The rotor is stationary.
- Press the *[Lid]* button.
    - ➔ The lid unlocks by means of a motor.
    - The light on the *[Lid]* button goes out.
    - The 'Lid open' indicator appears.

**Closing the lid**



**CAUTION**

**Danger of crushing when closing the lid.**

Danger of crushing fingers when the closing motor pulls the lid against the seal.

- No parts of the body should be in the hazard zone of the lid when closing the lid.
- To close the lid, press on the lid from above.



**NOTICE**

**Damage to the device caused by the lid slamming.**

- Close the lid slowly.
- Do not slam the lid.

- Personnel:**
- Trained user

- Close the lid and press down gently on the handle strip.
  - The lid locks using a motor.
  - The [Lid] button lights up.
  - The 'Lid closed' indicator appears.

## 6.2 Removing and installing the rotor

### Removing the rotor with a clamping nut




Fig. 28: Rotor installation and removal

- 1 Groove
- 2 Motor shaft
- 3 Driver
- 4 Hole

#### Personnel:

- Trained user
1. → Open the lid.
  2. → Loosen the rotor clamping nut using the supplied spanner.
    - After passing the working point for lifting the rotor, the rotor detaches from the cone of the motor shaft (2).
  3. → Turn the clamping nut until the rotor can be lifted off the motor shaft.
  4. → Remove the rotor.

### Installing the rotor with a clamping nut

#### Personnel:

- Trained user
- The lid is open.
1. → Clean the motor shaft (2) and rotor hole (4).
  2. → Lightly grease the motor shaft (2), see ➔ Chapter 8.2 'Cleaning and disinfection instructions' on page 42.
  3. → Place the rotor vertically on the motor shaft (2).
    - The driver (3) of the motor shaft must be in the groove (1) of the rotor. The orientation of the groove is marked on the rotor.
  4. → Hand-tighten the rotor clamping nut using the supplied spanner.
  5. → Check that the rotor is firmly seated.

## 6.3 Inserting and removing buckets

### Inserting buckets



#### NOTICE


Damage to the device due to imbalances caused by incorrect loading of the rotor.

- Load all swing-out rotor locations with the same buckets.



Buckets marked with the number of the rotor location may only be used there.

Buckets marked with a set number may only be used together.




1. Check that the rotor is firmly seated.
2. Grease the trunnions (3).
3. Insert the bucket (1) into the rotor from above. The trunnions (3) must be in the grooves (2).
4. Push the bucket (1) down as far as it will go.

**Removing the bucket**

- Pull the bucket (1) vertically upwards out of the rotor.

## 6.4 Inserting and removing adapters

**Inserting**



**the adapter**

- Insert the adapter (1) vertically into the bucket (2) from above.

**removing**

- Remove the adapter (1) vertically upwards out of the bucket (2).

## 6.5 Loading

**Filling centrifuge tubes**



**WARNING**

**Risk of injury from contaminated sample material.**

Contaminated sample material escapes from the sample tube during centrifugation.

- Use centrifuge tubes with special screw caps for hazardous substances.
- For risk group 3 and 4 materials, use a biosafety system in addition to the sealable centrifuge tubes (see WHO's 'Laboratory Biosafety Manual').



**NOTICE**

**Damage to the device due to highly corrosive substances.**

Highly corrosive substances may impair the mechanical strength of rotors, buckets and accessories.

- Do not centrifuge highly corrosive substances.



*Standard glass centrifuge tubes can be loaded up to RCF 4000 (DIN 58970 part 2).*

**Personnel:**

- Trained user





➔ Fill centrifuge tubes outside the centrifuge.

The maximum capacity of the centrifuge tubes specified by the manufacturer must not be exceeded.

With angle rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

It must be ensured that there is a uniform fill level in the tubes in order to keep the weight differences in the centrifuge tubes as low as possible.

**Loading the swing-out rotors**



**The following must be observed when using blood bags:**

**Personnel:**

- Trained user

1. ➔ Check that the rotor is firmly seated.

2. ➔ The centrifuge tubes must be distributed symmetrically across all rotor locations.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

No liquid must be allowed to enter the buckets and the centrifuging chamber when loading the buckets and swinging them out during the centrifugation run.

For containers with rubber inserts, there must always be the same number of rubber inserts under the centrifuge tubes.


All rotor locations must be filled with the same buckets. Certain buckets are marked with the number of the rotor location. The buckets must only be inserted in the corresponding rotor location.

Buckets marked with a set number (for example S001/4) must only be used in the set.

1. ➔ The differences can be compensated with balancing weights if the buckets are not filled with the same weight.

2. ➤ Empty buckets can be loaded with compensating inserts if there are insufficient blood bag systems available to load the rotor fully.
3. ➤ If necessary, fine-tune with the tare weights supplied.

**Loading the angle rotors**



**Personnel:**

- Trained user

1. ➤ Check that the rotor is firmly seated.
2. ➤ The centrifuge tubes must be distributed evenly over all locations on the rotor.

No liquid must be allowed to enter the rotor and the centrifuging chamber when loading the rotor.

With rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

**6.6 Opening and closing the biosafety system**

**6.6.1 Explanation**

The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.

Centrifuge tubes with special screw caps for hazardous substances must always be used.

For materials of risk group 3 and 4, a biosafety system must be used in addition to the sealable centrifuge tubes (see the World Health Organisation's "Laboratory Biosafety Manual").

In a biosafety system, a bioseal (sealing ring) prevents droplets and aerosols from escaping.

If the bucket of a biosafety system is used without the lid, the sealing ring must be removed from the bucket to prevent damage to the sealing ring during the centrifugation run.

Damaged biosafety systems are no longer microbiologically tight.

If no biosafety system is used, a centrifuge is not microbiologically tight in the sense of the EN / IEC 61010-2-020 standard.

**Storage of biosafety systems**

Biosafety systems must only be stored with the lid open to avoid damage to the sealing rings during storage.

## 6.6.2 Lid with spring-type lock




Fig. 29: Biosafety system

### Closing

1. ➤ Put on the lid.
2. ➤ Fold down both clamping brackets until they are under the lugs of the bucket.

### Opening

1. ➤ Fold up both clamping brackets until they are above the lugs of the bucket.
2. ➤ Remove the lid from the rotor.

## 6.7 Packing instructions, HettLiner


### Packing before centrifugation



Ensure that the plastic insert cannot tip over when loading and unloading the inserts (if necessary, use loading aid 4509).



1. ➤ Insert the blood bag (1) into the insert (2).



- 2.** → Hold the blood bag by the connectors (1) and push the support plate (2) on the outer side of the blood bag into the insert from top to bottom.

Ensure that the lower edge of the support plate stands on the floor as completely as possible.



- 3.** → Fold the support plate outwards and press it down until the folded edge of the support plate is level with the liquid level of the blood bag. The upper edge of the support plate must not project too far out of the insert during centrifugation due to the risk of its becoming caught in the rotor arms.

Note the position of the loop (1) so that it can be reached after centrifugation.



- 4.** → If available, fold empty satellite bag(s) and pack differently depending on the appropriate accessories and filling volume of the blood bag. It is advantageous to fold the satellite bags and pack them on the outside between the folded-down support plate and the outer wall of the insert.

Ensure that the silicone plate does not slip when doing this.




If necessary, the silicone plate can be retained and held against the loop when packing the satellite bag.

The position of the loop must be checked afterwards.



5.  Place the connections over the support plate so that the valves cannot break.  
Ensure that the hoses do not project out of the insert.  
Stow hose sections projecting over the edge of the insert between the folded-down support plate and the insert wall.
6.  If necessary, balance weights should be placed between the folded-down support plate and the receptacle wall.







#### Unpacking after centrifugation

1.  Pull the satellite bag out of the insert while holding the silicone plate in place with one hand.
2.  Slowly pull out the folded-down part of the support plate at the loop provided for this purpose.  
Return the support plate to its original shape in a controlled manner. The folded-down part of the support plate may spring back and mix blood components.
3.  Remove the remaining blood bag either together with the support plate or after removing the support plate from the insert.

## 6.8 Centrifugation



### 6.8.1 Centrifugation in continuous operation


#### Personnel:


- Trained user
1.  Press the *[t]* button repeatedly until the input field of the '*t/min:*' parameter is reverse highlighted.
  2.  Select the value 0 using the rotary knob.
  3.  Press the *[t]* button repeatedly until the input field of the '*t/sec:*' parameter is reverse highlighted.
  4.  Select the value 0 using the rotary knob.
    - The input field displays '*---:--*'.
  5.  Press the *[START]* button.
    - The centrifugation run is started.  
The '*Rotation*' indicator is lit up while the rotor is turning.  
The timing starts at 00:00.  
The rotor speed or the resulting RCF value, the temperature in the centrifuging chamber and the elapsed time are displayed during the centrifugation run.
  6.  Press the *[STOP]* button to cancel the centrifugation run.  
Ramp-down takes place with the set ramp-down parameters.
    - '*OPEN OEFFNEN*' is displayed.


### 6.8.2 Centrifugation with time preselection

#### Personnel:

- Trained user
1.  Press the *[t]* button repeatedly until the input field of the '*t/min:*' parameter is reverse highlighted.
  2.  Use the *[Rotary knob]* to set the desired value.

3.  Press the *[t]* button repeatedly until the input field of the '*t/:sec*' parameter is reverse highlighted.


4.  Use the *[Rotary knob]* to set the desired value.

5.  Press the *[START]* button.

➔ The centrifugation run is started.

The '*Rotation*' indicator is lit up while the rotor is turning.

The rotor speed or the resulting RCF value, the temperature in the centrifuging chamber and the remaining time are displayed during the centrifugation run.

6.  Ramp-down takes place with the selected ramp-down parameters after the time has elapsed or if the centrifugation run is cancelled by pressing the *[STOP]* button.


➔ '*OPEN OEFFNEN*' is displayed.

### 6.8.3 Changing settings during centrifugation

The runtime, speed, relative centrifugal force (RCF/RZB), ramp-up and ramp-down parameters, and the temperature (only for devices with cooling) can be changed during centrifugation.

These parameters can only be changed individually and one after the other.

1.  Change the value of the desired parameter using the *[Rotary knob]*

2.  Press the *[START]* button.

➔ The values of the current program are copied to program location '*---*' and updated with the changed value.

The original program is not overwritten.

## 6.9 Quick stop function

### Personnel:

- Trained user

 Press the *[STOP]* button twice.

➔ The '*STOP*' indicator flashes.

Ramp-down with brake level "R9" (shortest ramp-down time) is displayed and executed.

If brake level "R0" was selected, the ramp-down time is longer than with brake level "R9" for technical reasons.

## 7 Software operation

### 7.1 Key switches

The keys must be stored in such a way that they are protected from unauthorised access.

Key position	Function
Left key position	' <i>LOCK 1</i> ' is displayed. Programs can only be retrieved, not changed.
Right key position	' <i>LOCK 2</i> ' is displayed. No programs can be retrieved or changed.

Key position	Function
Middle key position	no status indicator. No program lock. Programs can be retrieved and changed.

## 7.2 Centrifugation parameters

### 7.2.1 Ramp-up and ramp-down parameters



The set ramp-up and ramp-down parameters are displayed.

x: 1-9 = ramp-up level, t = ramp-up time

y: R1-R9, B1-B9 = brake level, R0 = unbraked ramp-down, t = ramp-down time, n<sup>o</sup> = brake cut-off speed

#### Ramp-up level

1. ➤ Press the *[Ramp-up parameters]* button repeatedly until the 'Ramp-up level' parameter or the 'Ramp-up time' parameter is displayed and the input field is reverse highlighted.

2. ➤ Use the *[Rotary knob]* to set the desired level.

#### Ramp-up time

1. ➤ Press the *[Ramp-up parameters]* button repeatedly until the 'Ramp-up time min:sec' parameter is displayed and the input field is reverse highlighted.

2. ➤ Use the *[Rotary knob]* to set the desired level.

If a ramp-up time is set that is longer than the runtime, the centrifugation run ends before the set speed is reached.

#### Brake level

1. ➤ Press the *[Ramp-down parameters]* button repeatedly until the 'Ramp-down level' parameter or the 'Ramp-down time' parameter is displayed and the input field is reverse highlighted.

2. ➤ Use the *[Rotary knob]* to set the desired level.

B-brake levels can only be set for special rotors.

#### Ramp-down time

No ramp-down time can be set if a brake cut-off speed is set.

1. ➤ Press the *[Ramp-down parameters]* button repeatedly until the 'Ramp-down level min:sec' parameter is displayed and the input field is reverse highlighted.

2. ➤ Use the *[Rotary knob]* to set the desired level.

#### Brake cut-off speed

1. ➤ Press the *[Ramp-down parameters]* button repeatedly until the 'n<sup>o</sup>/RPM' parameter is displayed and the input field is reverse highlighted.

2. ➤ Use the *[Rotary knob]* to set the desired level.

### 7.2.2 Runtime



*For continuous operation, the minutes and seconds must be set to zero.*

*Continuous operation is indicated in the indicator by means of the '---:--' symbol.*



1. ➤ Press the *[t]* button repeatedly until the input field of the 't/min:.' parameter is reverse highlighted.

2. ➤ Use the *[Rotary knob]* to set the desired value.



3. ➤ Press the *[t]* button repeatedly until the input field of the 't/:sec' parameter is reverse highlighted.

4.  Use the *[Rotary knob]* to set the desired value.

### 7.2.3 Speed, RPM


1.  Press the *[n]* button repeatedly until the 'RPM' parameter is displayed and the input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the desired value.

### Maximum rotor speed indicator



1.  Press the *[n]* button repeatedly until the 'RPM' parameter is displayed and the input field is reverse highlighted.
2.  Press and hold the *[n]* button.
  - The maximum rotor speed (n-max rotor) is displayed.

### 7.2.4 Integral RCF

Integral RCF is a measure of the sedimentation effect ( $\int n^2 dt$ ). The value is used to compare centrifugation runs.

-  Press and hold the *[Integral RCF]* button.
  - 'Integral RCF' is displayed.

### 7.2.5 temperature

1.  Press the *[Temperature and centrifuging radius]* button repeatedly until the 'T/C°' parameter is displayed and the input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the desired value.

### 7.2.6 Relative centrifugal force, RCF

The relative centrifugal force RCF is dependent on the speed and the centrifuging radius.

The relative centrifugal force RCF is stated as a multiple of the acceleration due to gravity (g).

The relative centrifugal force RCF is a dimensionless numerical value and is used to compare the separation and sedimentation performance.

$$RCF = \left( \frac{RPM}{1000} \right)^2 * r * 1,118$$



$$RPM = \sqrt{\frac{RCF}{r * 1,118}} * 1000$$

RCF = Relative Centrifugal Force


RPM = speed


r = centrifuging radius in mm = distance from the centre of the axis of rotation to the bottom of the centrifuge tube.

### 7.2.7 Setting the relative centrifugal force (RCF/RZB)

1.  Press the *[RCF]* button repeatedly until the 'RCF/RZB' parameter is displayed and the input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the desired value.

### Maximum rotor RCF indicator

1.  Press the *[RCF]* button repeatedly until the 'RCF/RZB' parameter is displayed and the input field is reverse highlighted.

2.  Press and hold the *[RCF]* button.
  - ➔ The maximum RCF of the rotor (RCF-max rotor) is displayed.

## 7.2.8 Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm<sup>3</sup>

The density of the substances or mixtures of substances must not exceed 1.2 kg/dm<sup>3</sup> during centrifugation at maximum speed. The speed must be reduced for substances or substance mixtures with a higher density. The permissible speed can be calculated using the following formula:

$$\text{Reduzierte Drehzahl } (n_{red}) = \sqrt{\frac{1,2}{\text{höhere Dichte (kg/dm}^3)}} * \text{maximale Drehzahl (RPM)}$$

For example: Maximum speed 4000 RPM, density 1.6 kg/dm<sup>3</sup>

$$n_{red} = \sqrt{\frac{1,2(\text{kg/dm}^3)}{1,6(\text{kg/dm}^3)}} * 4000 \text{ RPM} = 3464 \text{ RPM}$$

If, in exceptional cases, the maximum load indicated on the bucket is exceeded, the speed must also be reduced. The permissible speed can be calculated using the following formula:



$$\text{Reduzierte Drehzahl } (n_{red}) = \sqrt{\frac{\text{maximale Beladung (g)}}{\text{tatsächliche Beladung (g)}}} * \text{maximale Drehzahl (RPM)}$$

For example: Maximum speed 4000 RPM, maximum load 300 g, actual load 350 g

$$n_{red} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} * 4000 \text{ RPM} = 3703 \text{ RPM}$$

Please contact the manufacturer if you are not sure.




## 7.2.9 Centrifuging radius

1.  Press the *[Temperature and centrifuging radius]* button repeatedly until the 'r/mm' parameter is displayed and the input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the desired value.

Changing the radius, automatically adjusts the RCF/RZB value: this is indicated by flashing.

## 7.3 Programming

### 7.3.1 Opening or loading programs

1.  Use the *[PROG]* button to select the 'PROG no.' parameter. The input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the desired program location.
3.  Press the *[RCL]* button.
  - ➔ The centrifugation data of the desired program location is displayed.

### 7.3.2 Entering or changing programs

1.  Set the desired parameters.

2. Use the *[PROG]* button to select the 'PROG no.' parameter. The input field is reverse highlighted.
3. Use the *[Rotary knob]* to set the desired program location.  
If the program location indicator flashes, this program location is already assigned centrifugation data. In this case, set a free program location, or overwrite the centrifugation data by continuing.
4. Press the *[STO]* button.
  - Settings are stored in the desired program location.
5. Press the *[STO]* button twice.
  - Centrifugation data that is already stored will be overwritten.

### 7.3.3 Automatic buffer

The buffer includes the program locations "----" and 90 to 99.

Changed centrifugation data is automatically saved to program location "----" every time a centrifugation run is started.

The changed centrifugation data, of the last 11 centrifugation runs, are stored in the buffer and can be retrieved.

## 7.4 Rotor detection

- Rotor detection is performed after starting a centrifugation run.
- If the rotor has been changed, the centrifugation run is cancelled after rotor detection. The rotor code (R) and the maximum rotor speed (n-max) of the newly detected rotor are displayed.
- If the maximum speed of the rotor used is less than the set speed, the speed is limited to the maximum rotor speed.

## 7.5 Cooling (for centrifuges with cooling)

### 7.5.1 Instructions, cooling

For centrifuges with the heating/cooling option, the temperature setpoint can be adjusted from -20 °C to +90 °C. If the actual temperature deviates from the set temperature by more than 5 °C, this is signalled by a flashing temperature value indicator.

The lowest achievable temperature is rotor dependent.

### 7.5.2 Standby cooling

The centrifuging chamber is cooled to the preselected temperature when the rotor is at a standstill and the lid is closed. The display shows the temperature setpoint.

### 7.5.3 Precooling the rotor

For rapid precooling of the unloaded rotor and accessories, we recommend a centrifugation run with the continuous operation settings and a speed of approx. 20% of the maximum rotor speed.

## 7.6 Heating (for centrifuges with heating)

During the centrifugation run, the centrifuging chamber is heated to the preselected temperature if required. The heating is switched off when the rotor is at a standstill.

**CAUTION****Danger of burns from hot surfaces.**

The surface temperature of the heating element in the centrifuging chamber can be up to 500 °C or 932 °F.

- Do not touch the heating element.
- The inside of the lid heats up during a centrifugation run with a very high temperature (e.g. +90 °C). In this case, do not touch the inside of the lid.

**NOTICE****Damage to plastic buckets due to excessive temperature**

- Plastic buckets may only be used at temperatures up to a maximum of 40 °C or 104 °F.










## 7.7 Machine Menu

### 7.7.1 Querying system information

The following system information can be queried:




- Centrifuge model
- Maximum speeds of the various rotor codes
- Centrifuge program version
- Frequency converter type
- Program version for the frequency inverter

The rotor is stationary.

1.  Press and hold the *[t]* button.
  - ➔ The audible signal 'SOUND / BELL' [sounds] after 8 seconds.
2.  Press the *[t]* button.
  - ➔ The operating hours 'CONTROL:' are displayed.
3.  Press the *[t]* button.
  - ➔ The date and time are displayed.
4.  Press the *[t]* button.
  - ➔ The machine and cooling version 'VERS 12 °C / \* 03' is displayed.
5.  Press the *[t]* button.
  - ➔ The operating hours of the frequency converter 'FC/CCI XX h' are displayed.
6.  Press the *[t]* button.
  - ➔ The type of frequency converter 'FU/CCI' is displayed.
7.  Press the *[t]* button.
  - ➔ The program version of the frequency converter 'FU/CCI - S.' is displayed.
8.  Press the *[t]* button.
  - ➔ The program version of the power supply circuit board '°C / \* - S. 01.07' is displayed.
9.  Press the *[STOP/OPEN]* button to exit the menu

## 7.7.2 Querying operating hours

The rotor is stationary.

1.  Open the lid.
2.  Press and hold the *[t]* button.
  - ➔ 'SOUND / BELL XXX' is displayed after 8 seconds.
3.  Press the *[t]* button.
  - ➔ 'CONTROL:' and the operating hours are displayed.

The operating hours indicator goes off automatically after 10 seconds.

## 7.7.3 Audible signal





### 7.7.3.1 General

The audible signal sounds after the following settings:

OFF	<ul style="list-style-type: none"> <li>■ after a problem occurs in the 2 s interval.</li> </ul>
ON1	<ul style="list-style-type: none"> <li>■ after a problem occurs in the 2 s interval.</li> <li>■ after completion of the centrifugation run and rotor standstill in the 30 s interval.</li> </ul>
ON2	<ul style="list-style-type: none"> <li>■ after a problem occurs in the 2 s interval.</li> <li>■ after completion of the centrifugation run and rotor standstill in the 30 s interval.</li> <li>■ every time the button is pressed.</li> </ul>




Opening the lid or pressing any button stops the audible signal.

### 7.7.3.2 Setting an audible signal


1.  Open the lid.
2.  Press and hold the *[t]* button.
  - ➔ 'SOUND / BELL ON1', 'SOUND / BELL ON2' or 'SOUND / BELL OFF' is displayed after 8 seconds.
3.  Use *[Rotary knob]* to set 'OFF', 'ON1' or 'ON2'.
4.  Press the *[START]* button.
  - ➔ The setting is stored.
  - \*\*\* OK \*\*\* is displayed briefly.

## 7.7.4 Centrifugation data displayed after switching on

The centrifugation data of program 1 or the last program used is displayed after switching on.







1.  Set the mains switch to *[I]*.
2.  Press the *[STOP]* button at the first visual change in the display (inverse display).
  - ➔ 'PROGRAM 1, LAST PROGRAM' is displayed.
3.  Use the *[Rotary knob]* to set the desired function.



4.  Press the *[START]* button.
  - The settings are stored.
  - \*\*\* OK \*\*\** is displayed briefly.

## 7.7.5 Setting the date and time

The rotor is stationary.

1.  Open the lid.
2.  Press and hold the *[t]* button.
  - *'SOUND / BELL'* is displayed after 8 seconds.
3.  Press the *[t]* button twice.
  - The date and time are displayed
    - a: Year
    - mon: Month
    - d: Day
    - h: Hours
    - min: Minutes
4.  Press the *[Temperature and centrifuging radius]* button repeatedly until the desired parameter is displayed and the input field is reverse highlighted.
5.  Use the *[Rotary knob]* to set the desired value.
6.  Press the *'Start'* button.
  - The settings are stored.
  - \*\*\* OK\*\*\** is displayed briefly.

## 7.8 Program links

### 7.8.1 Linking programs or changing a program link












*Program linking is only possible with programs where ramp-up and brake levels are set.*

*The programs must be saved in the desired order before linking, by means of either program entry or program retrieval.*





*The program locations must be consecutive (e.g. program locations 10+11+12).*

#### Linking programs






1.  Use the *[PROG]* button to select the *'PROG no.'* parameter. The input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the program location of the initial program (XX+).
3.  Press the *[RCL]* button.
  - The centrifugation data of the desired program location is displayed
4.  Press the *[PROG]* button twice.
  - Parameter PR-PART is selected.
  - The input field is reverse highlighted.


5.  Press the *[STO]* button twice.
  - The program is linked and the program number of the next program location (+XX+) is displayed.
6.  Press the *[RCL]* button twice.
  - The centrifugation data of the desired program location is displayed
7.  Press the *[STO]* button twice.
  - The program is linked and the program number of the next program location (+XX+) is displayed.
8.  Repeat the last two steps until all programs are linked.
9.  Press the *[PROG]* button.
  - The program number of the end program (+XX) is displayed.

### Changing a program link

1.  Open the desired program.
2.  Change the desired parameter.
3.  Save the changed centrifugation data to the same program location again.
  - Saving removes the program link.
4.  Link programs again.

### 7.8.2 Centrifugation run with program link

1.  Press the *[PROG]* button twice.
  - Parameter PR-PART is selected.  
The input field is reverse highlighted.
2.  Use the *[Rotary knob]* to set the program location of the initial program (XX+).
3.  Press the *[RCL]* button.
  - The centrifugation data of the desired program location is displayed
4.  Press the *[START]* button.
  - The centrifugation run is started.  
The 'Rotation' indicator appears as long as the rotor is turning.  
The ramp-up and brake levels of the program link are displayed.
    - Initial program (XX+)
      - x: Ramp-up level of the initial program x
    - Follow-up program (+XX+)
      - x: Ramp-up level of the follow-up program x
    - End program (+XX)
      - x: Ramp-up level of the end program
      - y: Brake level of the end program
5.  Ramp-down takes place with the brake level of the end program after the time in the end program has elapsed.  
If the centrifugation run is cancelled by pressing the *[STOP]* button, ramp-down takes place at the brake level of the program currently running.



### 7.8.3 Deleting program links

1. Use the *[PROG]* button to select the 'PROG no.' parameter. The input field is reverse highlighted.
2. Use the *[Rotary knob]* to set the program location of the initial program (XX+).
3. Press the *[RCL]* button.
  - The centrifugation data of the desired program location is displayed
4. Press the *[PROG]* button twice.
  - Parameter 'PR-PART' is displayed.  
The input field is reverse highlighted.
5. Press the *[STO]* button twice.
6. Press the *[PROG]* button.

## 8 Cleaning and care

### 8.1 Overview table

Chap.	Task to execute	if required	daily	weekly	Annually	Page
<b>8</b>	<b>Cleaning and care</b>					41
<b>8.3</b>	<b>Cleaning</b>					42
8.3	Cleaning the device		X			42
8.3	Cleaning the biosafety systems			X		42
8.3	Cleaning the accessories			X		43
<b>8.4</b>	<b>Disinfection</b>					43
8.4	Disinfecting the device	X				43
8.4	Disinfecting the accessories	X				43
<b>8.5</b>	<b>Maintenance</b>					44
8.5	Greasing the rubber seal of the centrifuging chamber			X		44
8.5	Greasing the rubber seal of the biosafety system			X		44
8.5	Trunnion greasing			X		44
8.5	Checking the accessories			X		44
8.5	Checking the biosafety system			X		44
8.5	Inspecting the centrifuging chamber for damage				X	44
8.5	Greasing the motor shaft				X	44

Chap.	Task to execute	if required	daily	weekly	Annually	Page
8.5	Accessories with a limited service life	X				45
8.5	Replacing centrifuge tubes	X				45

## 8.2 Cleaning and disinfection instructions



### DANGER

Risk of contamination for the user due to inadequate cleaning or failure to observe the cleaning instructions.

- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Observe laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.

- The device and its accessories must not be cleaned in dishwashers.
- Only perform hand cleaning and liquid disinfection.
- The water temperature must not exceed 25 °C.
- To prevent any corrosion due to use of detergents or disinfectants, it is essential to follow the special application instructions provided by the manufacturers of the detergent or disinfectant.

### Disinfectant:

- Surface disinfectant (not disinfectant for hands or instruments)
- Ethanol as the sole active substance.  
Do not use an ethanol-propanol mixture to disinfect the viewing window in the lid of the device.
- Concentration is not less than 30 %
- pH: 6 – 8
- Non-corrosive

## 8.3 Cleaning

### Cleaning the device

1. ➤ Open the lid.
2. ➤ Switch off the device and disconnect it from the power supply.
3. ➤ Remove accessories.
4. ➤ Clean the centrifuge housing and the centrifuging chamber with soap or a mild detergent and a damp cloth.
5. ➤ Remove any detergent residues with a damp cloth after using detergents.
6. ➤ The surfaces must be dried immediately after cleaning.
7. ➤ Dry the centrifuging chamber with an absorbent cloth if condensation forms.

### Cleaning the biosafety systems

1. ➤ Clean the biosafety system using the detergent and a damp cloth.

2. ➤ Remove any detergent residues with a damp cloth after using detergents.
3. ➤ Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

### Cleaning the accessories

1. ➤ Clean the accessories using the detergent and a damp cloth.
2. ➤ Remove any detergent residues with a damp cloth after using detergents.
3. ➤ Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

## 8.4 Disinfection



*Disinfection must always be preceded by cleaning of the components concerned.*

*See → Chapter 8 'Cleaning and care' on page 41*



*Disinfectant concentration and application time according to the manufacturer's instructions.*

### Disinfecting the device



#### CAUTION

**Risk of injury due to ingress of water or other liquids.**

- Protect the device against external liquids.
- Do not disinfect the device using spray.

1. ➤ Open the lid.
2. ➤ Switch off the device and disconnect it from the power supply.
3. ➤ Remove accessories.
4. ➤ Clean the housing and centrifuging chamber using disinfectant.
5. ➤ Remove any disinfectant residues with a damp cloth after using disinfectants.
6. ➤ The surfaces must be dried immediately after cleaning.

### Disinfecting the accessories

1. ➤ Disinfect the accessories using the disinfectant.
2. ➤ Wet all cavities with bubble-free disinfectant.
3. ➤ Remove the disinfectant residues or leave them to dry after using disinfectants.

### Autoclaving

The following accessories may be autoclaved at 121 °C / 250 °F (20 min):

- Swing-out rotors
- Aluminium angle rotors
- Metal buckets
- Lid with bioseal
- Inserting

No statement can be made about the resulting degree of sterility.

The lids of the rotors and bucket must be removed before autoclaving.

Autoclaving accelerates the ageing of materials. It may cause changes to colours. After autoclaving, the rotors and accessories are to be visually inspected for damage and any damaged parts are to be replaced immediately.

The sealing ring in question is to be replaced if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

The sealing rings must be replaced after autoclaving to ensure the tightness of the biosafety systems.

## 8.5 Maintenance

### Greasing the rubber seal of the centrifuging chamber

→ Rub the sealing ring lightly with a rubber care product.

### Greasing the rubber seal of the biosafety system

→ Rub the sealing ring lightly with a rubber care product.

### Trunnion greasing

1. → Remove accessories.
2. → Clean the trunnions.
3. → Remove any detergent residues with a damp cloth after using detergents.
4. → Grease the trunnions and suspension with Hettich Tubenfett 4051.
5. → Excess grease in the centrifuging chamber must be removed.

### Checking the accessories

1. → The accessories are to be checked for wear and corrosion damage.
2. → Check that the rotor is firmly seated.

### Checking the biosafety system

1. → Visually check all parts of the biosafety system for damage.
2. → Check the correct installation position of the sealing ring(s) of the biosafety system.
3. → Replace the damaged parts of the biosafety system.
4. → Replace the sealing ring in question immediately if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

### Inspecting the centrifuging chamber for damage

→ Check the centrifuging chamber for damage.

### Greasing the motor shaft

1. → Remove accessories.
2. → Clean the motor shaft.
3. → Remove any detergent residues with a damp cloth after using detergents.
4. → Grease the motor shaft with Hettich Tubenfett 4051.
5. → Excess grease in the centrifuging chamber must be removed.

**Accessories with a limited service life**

The use of certain accessories is time-limited. For safety reasons, the accessories must no longer be used when either the maximum number of permissible run cycles marked on them or the expiry date marked on them has been reached.

- The maximum permissible number of run cycles or the expiry date can be seen marked on the accessories.
- The centrifuge is equipped with a cycle counter.

**Replacing centrifuge tubes**



**CAUTION**

**Risk of injury from broken glass.**

Broken glass may cause glass splinters and contaminated liquids to be found inside the centrifuge.

- Wear cut-resistant gloves.
- Wear protective goggles and a face mask.

Broken parts of the tube, glass splinters and spilled centrifuge material must be removed completely in the event of leakage or if a centrifuge tube breaks. Glass splinters that are not removed will cause further glass breakage.

The rubber inserts and the plastic sleeves of the rotors must be replaced after a glass breakage.

Disinfection must be carried out if the material is infectious.


## 9 Troubleshooting

### 9.1 Fault description




Customer service must be notified if the fault cannot be rectified based on the fault table. State the centrifuge type and serial number. Both numbers can be seen on the type plate of the centrifuge.

\* Error number does not appear on the display.

Fault description	Cause	Remedy
no display	No power. Triggering of the overcurrent protection fuse (for centrifuges 5005-08 only).	<ul style="list-style-type: none"> <li>■ Check the supply voltage.</li> <li>■ The mains switch is in switch position <i>///</i></li> </ul>
TACHO - ERROR 01, 02	Tacho defective. Motor, inverter, electronics defective.	<ul style="list-style-type: none"> <li>■ Open the lid.</li> <li>■ Set the mains switch to <i>[0]</i>.</li> <li>■ Wait at least 10 seconds.</li> <li>■ Turn the rotor vigorously by hand.</li> <li>■ Set the mains switch to <i>///</i>. The rotor must rotate while switching on.</li> </ul>
IMBALANCE / UNWUCHT	The rotor is unevenly loaded.	<ul style="list-style-type: none"> <li>■ Open the lid or hatch.</li> <li>■ Check the loading of the rotor.</li> <li>■ Repeat the centrifugation run.</li> </ul>
CONTROL - ERROR 04, 06-09	Lid lock error.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
N > MAX 05	Overspeed error	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>

Fault description	Cause	Remedy
N < MIN 13	Underspeed error.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
ROTORCODE 10	Rotor coding error.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
MAINS INTERRUPT	Loss of mains power during the centrifugation run. The centrifugation run was not completed.	<ul style="list-style-type: none"> <li>■ Open the lid.</li> <li>■ Press the <i>[START]</i> button.</li> <li>■ If required: Repeat the centrifugation run.</li> </ul>
VERSION-ERROR 12	No conformity of the electronic components, error/defect in electronics.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
SER I/O - ERROR 30-38	Error/defect in interface.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
° C * - ERROR 50-56, 58	Error/defect in cooling.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
LOCK - ERROR 57	Error/defect in program lock.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
FU/CCI-ERROR 60-83	Error/defect in motor control.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
CONTROL - ERROR 26, 90-95, 97 - 99	Error/defect in control section.	<ul style="list-style-type: none"> <li>■ Perform a MAINS RESET.</li> </ul>
N > ROTOR MAX 96	Speed in the selected program greater than the maximum rotor speed.	<ul style="list-style-type: none"> <li>■ Check and correct the speed.</li> </ul>
	The rotor has been changed. The built-in rotor has a higher maximum speed than the previously used rotor. The rotor has not yet been recognised by the rotor detection.	<ul style="list-style-type: none"> <li>■ Set a speed up to the maximum speed of the previously used rotor. Press the <i>[START]</i> button to perform rotor detection.</li> </ul>
 The entire display lights up.	-	<ul style="list-style-type: none"> <li>■ Notify customer service.</li> </ul>

## 9.2 Perform a MAINS RESET

1.  Set the mains switch to *[0]*.
2.  Wait 10 seconds.
3.  Set the mains switch to *[I]*.

## 9.3 Emergency release

The lid cannot be unlocked by the motor in the event of a power failure. Emergency unlocking by hand must be performed.



### WARNING

Risk of electric shock due to maintenance and servicing work on live device.

- Disconnect the device from the mains before carrying out repairs and maintenance.





**WARNING**

**Danger of cutting and crushing due to moving rotor.**

- Do not open the lid until the rotor has stopped.

**Personnel:**

- Trained user

1. ▶ Look through the window in the lid to ensure that the rotor is stationary.
2. ▶ Insert the hex key horizontally into the hole (1) and turn clockwise half a turn until the lid opens.
3. ▶ Remove the hex key from the hole (1).
4. ▶ When the power is restored, press the [Lid] button so that the motorised cover lock returns to the home position (open).




Fig. 30: Emergency release

1 Hole

## 10 Disposal

### 10.1 General instructions



**The device can be disposed of via the manufacturer.**

A Return Material Authorisation (RMA) form must always be requested for a return.

If necessary, contact the Technical Service Department of the manufacturer.

- **Andreas Hettich GmbH & Co. KG**
- Föhrenstrasse 12
- 78532 Tuttlingen, Germany
- Phone: +49 7461 705 1400
- E-mail: [service@hettichlab.com](mailto:service@hettichlab.com)



**! WARNING**

**Risk of pollution and contamination for people and the environment.**

When disposing of the centrifuge, people and the environment may be polluted or contaminated by incorrect or improper disposal.

- Removal and disposal may be carried out only by a trained and authorised service personnel.

The device is intended for the commercial sector ("Business to Business" - B2B).

According to Directive 2012/19/EU, the devices may no longer be disposed of with household waste.

The devices are assigned to the following groups according to the Stiftung Elektro-Altgeräte Register (EAR (German foundation under civil law)):

- Group 1 (heat exchangers)



The crossed-out wheelie bin symbol indicates that the device must not be disposed of with household waste. Regulations governing disposal of such devices may differ in individual countries. If necessary, contact the supplier.



*Fig. 31: Household waste ban*

## 11 Index

### A

Akustisches Signal	
aktivieren/deaktivieren. . . . .	38
Allgemeine Sicherheitshinweise. . . . .	8
An- und Auslaufparameter. . . . .	33
Anschließen der Zentrifuge. . . . .	22
Aufstellen der Zentrifuge. . . . .	22
Auspacken. . . . .	21
Ausschalten. . . . .	24
Autoklavieren. . . . .	43

### B

Befüllen. . . . .	26
Beladen. . . . .	26
Betriebsstunden	
abfragen. . . . .	38
Bio-Sicherheitssystem	
prüfen. . . . .	44
reinigen. . . . .	42

### D

Datums- und Uhreinstellung. . . . .	39
Dauerlauf. . . . .	31
Desinfektion. . . . .	43
Drehzahl RPM. . . . .	34

### E

Einschalten. . . . .	23
Einstellung während Zentrifugation. . . . .	32
Entsorgung. . . . .	47
Ersatzteile. . . . .	19

### F

Fehlermeldungen. . . . .	45
--------------------------	----

### G

Gerät	
desinfizieren. . . . .	43
reinigen. . . . .	42
Gummidichtung	
fetten. . . . .	44

### I

Integrale Zentrifugalbeschleunigung	
Integral RCF. . . . .	34

### L

Lid	
the lid. . . . .	24
Lieferumfang. . . . .	19

### M

Motorwelle	
fetten. . . . .	44

### N

NETZ-RESET. . . . .	46
Nicht vorgesehene Zweckbestimmung. . . . .	7

### O

Originalersatzteile. . . . .	19
------------------------------	----

### P

Personalunterweisung. . . . .	8
Persönliche Schutzausrüstung. . . . .	7
Personnel qualifications. . . . .	7
Pflege	
Intervalle. . . . .	41
Programm	
ändern. . . . .	35
aufrufen. . . . .	35
eingeben. . . . .	35
laden. . . . .	35
Programmverknüpfung	
ändern. . . . .	40
erstellen. . . . .	39
löschen. . . . .	41
Zentrifugationslauf. . . . .	40

### R

Reinigung. . . . .	42
Reinigung und Desinfektion	
Hinweise. . . . .	42
Relative Zentrifugalbeschleunigung	
RCF. . . . .	34
Relative Zentrifugalbeschleunigung (RCF/RZB). . . . .	34
Rotor	
ausbauen. . . . .	25
beladen. . . . .	27, 28
einbauen. . . . .	25
Rotorerkennung. . . . .	36
Rücksendung. . . . .	19
Runtime. . . . .	33

### S

Schilder	
am Gerät. . . . .	14
auf der Verpackung. . . . .	14
Schleuderraum	
prüfen. . . . .	44
Schlüsselschalter. . . . .	32
Schutzausrüstung. . . . .	7
Sicherheitshinweise. . . . .	8
Storage conditions. . . . .	20
Störungsbehebung. . . . .	45
Symbole. . . . .	6
Systeminformationen	
abfragen. . . . .	37

### T

Tragzapfen	
fetten. . . . .	44
Transport condition. . . . .	20
Trouble shooting. . . . .	45
Typenschild. . . . .	13

**V**

Verantwortung des Betreibers. . . . .	8
Vorgesehene Zweckbestimmung. . . . .	6
Vorhersehbare Fehlanwendung. . . . .	7

**W**

Wartung. . . . .	44
Intervalle. . . . .	41

**Z****Zentrifugation**

im Dauerlauf. . . . .	31
mit höherer Stoffdichte. . . . .	35
mit Zeitvorwahl. . . . .	31
Zentrifugationsdaten nach Einschalten. . . . .	38

**Zentrifugiergefäße**

tauschen. . . . .	45
-------------------	----
















**Zubehör. . . . . 19**



















desinfizieren. . . . .	43
mit begrenzter Verwendungsdauer. . . . .	45
prüfen. . . . .	44
reinigen. . . . .	43

**Zwischenspeicher**

automatisch. . . . .	36
----------------------	----















# Rotoren und Zubehör / Rotors and accessories















4174	5052 + 5057								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>      max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g  mit Bioabdichtung / with bio-containment 9)								
	4213		4213-93			4214			
									
	---	---	---	---	---	---	---	---	
	 <sup>2)</sup>	 <sup>2)</sup>		 <sup>2)</sup>				---	
Kapazität / capacity	ml	6	7	4,5 - 5	5	2,7 - 3	7,5 - 8,2	5 - 10	---
Maße / dimensions $\varnothing$ x L	mm	12 x 82	12 x 100	11 x 92	12 x 75	11 x 66	15 x 92	16 x 100	---
Anzahl p. Rotor / number p. rotor		192	192	192	192	192	120	120	---
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	5003							
Radius / radius	mm	221							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	-1							

4174	5052 + 5057								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>      max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g  mit Bioabdichtung / with bio-containment 9)								
	4214		4214-93			4215		4216	4218
									
	---	---	---	---	---	---	---	---	
		 <sup>2)</sup>	 <sup>2)</sup>			 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>	
Kapazität / capacity	ml	10	15	4 - 5,5	4 - 7	25	30	50	100
Maße / dimensions $\varnothing$ x L	mm	15 x 102	17 x 100	15 x 75	16 x 75	24 x 100	25 x 110	34 x 100	44 x 100
Anzahl p. Rotor / number p. rotor		120	120	120	120	44	44	24	16
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	5003							
Radius / radius	mm	221							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	-1							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.



- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".




4174		5052 + 5057								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>     90°  max. Laufzyklen / max. cycles: 30000 (4500RPM) 60000 (4000 RPM) 80000 (3500 RPM)		  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g  mit Bioabdichtung / with bio-containment <sup>9)</sup>								
		4220		4222		4222-93		---		
								---		
		---	---	---	---	---	---	---	---	
								---		
Kapazität / capacity	ml	9 - 10	12	7	4 - 7	5	2,6 - 3,4	1 - 5	---	
Maße / dimensions	∅ x L	mm	16 x 92	17 x 100	12 x 100	13 x 100	12 x 75	13 x 65	13 x 75	---
Anzahl p. Rotor / number p. rotor		64	64	120	120	120	120	120	---	
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5003								
Radius / radius	mm	221								
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-1								

4174		5052 + 5057									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>     90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g  mit Bioabdichtung / with bio-containment <sup>9)</sup>									
		---		4223		---		4224		4225	
		---				---					
		---		---	---	---	---	---	---	---	---
		---				---					
Kapazität / capacity	ml	---	9	8	12	---	4	1,5	2		
Maße / dimensions	∅ x L	mm	---	14 x 100	16 x 125	16 x 101	---	10 x 88	11 x 38	11 x 38	
Anzahl p. Rotor / number p. rotor		---	100	100	100	---	252	336	336		
Drehzahl / speed	RPM	---	4500	4500	4500	---	4500	4500	4500		
RZB / RCF	<sup>2)</sup>	---	5003	5003	5003	---	5003	5094	5094		
Radius / radius	mm	---	221	221	221	---	221	225	225		
 9 (97%)	sec	125									
 9	sec	197									
Temperatur / temperature	°C <sup>1)</sup>	-1									

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
















4174	5052 + 5057							
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)								
	max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM)							
	max. Beladung / max. load: 1100 g							
	mit Bioabdichtung / with bio-containment <sup>9)</sup>							
	4226	4232	4249	SK 09.00	4238	4241	4245-A	
	Microtainer	---	---	---	5127	---	---	
					8)	8)		
	Kapazität / capacity ml	0,8	15	50	14	250	290	25
Maße / dimensions Ø x L mm	8 x 45	17 x 120	29 x 115	16,5 x 106	62 x 122	62 x 137	25 x 90	29 x 115
Anzahl p. Rotor / number p. rotor	312	92	24	92	4	4	32	32
Drehzahl / speed RPM	4500							
RZB / RCF <sup>2)</sup>	5003	5184	5094	5117	5003	5003	5184	
Radius / radius mm	221	229	225	226	221	221	221	269
9 (97%)	125							
9	197							
Temperatur / temperature °C <sup>1)</sup>	-1							


















4174	5052 + 5057				5052	---			
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)									
	max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM)						---	---	---
	max. Beladung / max. load: 1100 g						---	---	---
	mit Bioabdichtung / with bio-containment <sup>9)</sup>					---	---	---	
	4258			4258+4449		6322	---	---	---
							---	---	---
	0512	0554	4234-A	Corning	Corning	---	---	---	
						---	---	---	
		8)	8)	8)					
	Kapazität / capacity ml	750	650	750	500	250	---		
Maße / dimensions Ø x L mm	97 x 152	97 x 139	96 x 135	96 x 147	60 x 162	---			
Anzahl p. Rotor / number p. rotor	4				6	---			
Drehzahl / speed RPM	4500					---			
RZB / RCF <sup>2)</sup>	5184			5184	5003	---			
Radius / radius mm	229			229	221	---			
9 (97%)	125					---			
9	197					---			
Temperatur / temperature °C <sup>1)</sup>	-1					---			

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 3) Gefäß nur belastbar bis RZB 700
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 3) tube will not stand RCF values exceeding 700
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".































4174	4522-A								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g								
	4213		4213-93			4214			
									
	---	---	---	---	---	---	---	---	
	 2)	 2)		 2)					
Kapazität / capacity	ml	6	7	4,5 - 5	5	2,7 - 3	7,5 - 8,2	5 - 10	10
Maße / dimensions $\varnothing \times L$	mm	12 x 82	12 x 100	11 x 92	12 x 75	11 x 66	15 x 92	16 x 100	115 x 102
Anzahl p. Rotor / number p. rotor		192	192	192	192	192	120	120	120
Drehzahl / speed	RPM	4500							
RZB / RCF		5003							
Radius / radius	mm	221							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	-1							

4174	4522-A								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g								
	4214		4214-93		4215		4216	4218	
									
	---	---	---	---	---	---	---	---	
		 2)			 2)		 2)	 2)	
Kapazität / capacity	ml	15	15	4 - 5,5	4 - 7	25	30	50	100
Maße / dimensions $\varnothing \times L$	mm	17 x 100	17 x 100	15 x 75	16 x 75	24 x 100	25 x 110	34 x 100	44 x 100
Anzahl p. Rotor / number p. rotor		120		120		44		24	16
Drehzahl / speed	RPM	4500							
RZB / RCF		5003							
Radius / radius	mm	221							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	-1							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.




















- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.

4174	4522-A									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\sphericalangle$ 90  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g									
	4220		4222		4222-93		---			
							---			
	---	---	---	---	---	---	---	---		
								---		
Kapazität / capacity	ml	9 - 10	12	7	4 - 7	5	2,6 - 2,9	1 - 5	---	
Maße / dimensions	Ø x L	mm	16 x 92	17 x 100	12 x 100	13 x 100	12 x 75	13 x 65	13 x 75	---
Anzahl p. Rotor / number p. rotor		64		120		120		---		
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5003								
Radius / radius	mm	221								
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-1								

4174	4522-A									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM) max. Beladung / max. load: 1100 g									
	---	4223			---	4224		4225		
	---				---					
	---	---	---	---	---	---	---	---	---	
	---				---					
Kapazität / capacity	ml	---	9	8	12	---	4	1,5	2	
Maße / dimensions	Ø x L	mm	---	14 x 100	16 x 125	16 x 101	---	10 x 88	11 x 38	11 x 38
Anzahl p. Rotor / number p. rotor		---	100			---	252	336		
Drehzahl / speed	RPM	---	4500	4500	4500	---	4500	4500	4500	
RZB / RCF	<sup>2)</sup>	---	5003	5003	5003	---	5003	5094	5094	
Radius / radius	mm	---	221	221	221	---	221	225	225	
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-1								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.

4174		4522-A							
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)									
		max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM)							
		max. Beladung / max. load: 1100 g ---							
		---							
		4226	4232	SK 09.00	4249	4238	4241	4245-A	
									
		Microtainer	---	---	---	5127	---	---	---
						 8)	 8)		
Kapazität / capacity	ml	0,8	15	14	50	250	290	25	50
Maße / dimensions $\varnothing \times L$	mm	8 x 45	17 x 120	16,5 x 106	29 x 115	62 x 122	62 x 137	25 x 90	29 x 115
Anzahl p. Rotor / number p. rotor		312	92	138	24	4	32	32	
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	5003	5184	5117	5094	5003	5003	5003	5184
Radius / radius	mm	221	229	226	225	221	221	221	229
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	-1							

1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur



















2) Angaben des Röhrchenherstellers beachten.

8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature

2) Observe the tube manufacturer's instructions.

8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4174	4522-A					4524-A				
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>      max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)										
	max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM)					max. Laufzyklen / max. cycles: 15000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)				
	max. Beladung / max. load: 1100 g					max. Beladung / max. load: 2160 g				
	---					---				
	4258		4258+4449		6322	4529-AO <sup>10)</sup>	4529-AM <sup>10)</sup>	4529-AU <sup>10)</sup>		
										
0512		0554	4234-A	Corning	Corning	---	---	---		
 8)		 8)	 8)			 13)				
Kapazität / capacity	ml		750		650		750		500	
Maße / dimensions	Ø x L		mm		97 x 152		97 x 139		96 x 135	
Anzahl p. Rotor / number p. rotor	4		4		4		4		8	
Drehzahl / speed	RPM		4500		4500		4500		4500	
RZB / RCF	2)		5184		5184		5003		5683	
Radius / radius	mm		229		229		221		251	
 9 (97%)	sec		125		125		125		125	
 9	sec		197		197		197		197	
Temperatur / temperature	°C <sup>1)</sup>		-1-		-1-		-1-		10	

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 10) Darf nur in Gehänge 4524-A und beidseitig beladen verwendet werden. Bei Verwendung der Haltestifte darf die max. RZB 1000 nicht überschritten werden.




















	Position der Haltestifte	Abstand der Haltestifte vom Einsatzboden (innen)
4529-AO	oben	199,5 mm
4529-AM	Mitte	182,0 mm
4529-AU	unten	164,5 mm















- 13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 10) May only be used in bucket 4524-A and with both sides loaded. If the holding pins are used the max. RCF of 1000 may not be exceeded.

	Position of the holding pins	Distance of the holding pins from the bottom of the inserts (inside)
4529-AO	upper section	199.5 mm
4529-AM	middle section	182.0 mm
4529-AU	lower section	164.5 mm

















- 13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".













4174	4524-A		4546-A		4591-A		4595-C					
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  $\angle 90^\circ$ max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)							 max. 3500 RPM zulässig					
	max. Laufzyklen / max. cycles: 15000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)		max. Laufzyklen / max. cycles: 30000 (4500 RPM) 35000 (4000 RPM) 40000 (3500 RPM)		max. Laufzyklen / max. cycles: 15000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)		max. Laufzyklen / max. cycles: 4500 (3500 RPM) 15000 (3000 RPM)					
	max. Beladung / max. load: 2160 g		max. Beladung / max. load: 2200 g		max. Beladung / max. load: 2160 g		max. Beladung / max. load: 2570 g					
	4592-B		4559-A		4598-A	4592-B		4596-A				
												
4-fach / 4-times		1-fach / 1-times		4587-A	4-fach / 4-times	4584-A	3-fach / 3-times	4-fach / 4-times	1-fach / 1-times	4587-A	---	
 13)		 12)		 13)	 12)	 13)	 13)	 12)	---			
Kapazität / capacity	ml		500	750	---	450	---	450	500	750	---	max. 1600
Maße / dimensions $\varnothing \times L$	mm		---	---	---	---	---	---	---	---	---	---
Anzahl p. Rotor / number p. rotor	8		8		8		8		8		4	
Drehzahl / speed	RPM		4500									3500
RZB / RCF	2)		5683		5479		5705		5705		3328	
Radius / radius	mm		251		242		252		252		243	
 9 (97%)	sec		125									95
 9	sec		197									131
Temperatur / temperature	$^\circ\text{C}^{1)}$		10									-12

4174	4579-A	4579-A + 4255 / 4255-P <sup>4)</sup>									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  $\angle 90^\circ$ max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)											
	---	max. Laufzyklen / max. cycles: 30000									
	---	max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RPM)									
	---	---	---	---	4449	4430			---		
---	---	---	---					---			
---	4255 / 4255-P 4)	0512	4239	Corning		Nalgene®	Nunc®	---			
---	 11)	 5) 8)	 5) 8)	 5)	 5)	 5)	 5)	---			
Kapazität / capacity	ml		1000	750	1000	500	250	175	200	---	
Maße / dimensions $\varnothing \times L$	mm		98 x 138	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 139,2	60 x 130	---	
Anzahl p. Rotor / number p. rotor	4		4							---	
Drehzahl / speed	RPM		4500								
RZB / RCF	2)		5501								
Radius / radius	mm		243								
 9 (97%)	sec		125								
 9	sec		197								
Temperatur / temperature	$^\circ\text{C}^{1)}$		-11								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrencherstellers beachten.
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 5) 4255 nicht mit Deckel verschließbar
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 11) Maximale Beladung 800g. Bei einer Beladung über 800g muss die Drehzahl reduziert werden, siehe Beschriftung auf dem Becher. Berechnung der reduzierten Drehzahl siehe Kapitel "Zentrifugation von Stoffen oder Stoffgemischen mit einer höheren Dichte als 1,2 kg/dm<sup>3</sup>".
- 12) Ausgleichseinsatz. Handhabung siehe in Kapitel "Beladen des Rotors".
- 13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 11) Maximum load 800g. With a load higher than 800g the speed has to be reduced, see label on the bucket. Calculation of the reduced speed see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>".
- 12) Compensation insert. For handling see chapter "Loading the rotor".
- 13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".

4174	4579-A + 4255 / 4255-P <sup>4)</sup>									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800g (4500 RPM) 1000g (4020 RPM) 1200g (3670 RPM)									
	4432			4433			4434			
										
	---	---	---	---	---	---	---	---	---	---
			 2)	 2)			 2)	 2)	 2)	
Kapazität / capacity	ml	1,5	2,0	5	7	2,7 - 3	4,5 - 5	9	15	15
Maße / dimensions	∅ x L	mm	11 x 38	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	17 x 100	17 x 100
Anzahl p. Rotor / number p. rotor		168			120			76		
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5003			4935			5094		
Radius / radius	mm	221			218			225		
 9 (97%)	sec							125		
 9	sec							197		
Temperatur / temperature	°C <sup>1)</sup>	-11								

4174	4579-A + 4255 / 4255-P <sup>4)</sup>									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RPM)									
	---									
	4434									
										
	---	---	---	---	---	---	---	---	---	---
 5)		---						---	---	
Kapazität / capacity	ml	10	10	---	4 - 5,5	9 - 10	4 - 7	5 - 10	---	
Maße / dimensions	∅ x L	mm	15 x 102	16 x 80	---	15 x 75	16 x 92	16 x 75	16 x 100	---
Anzahl p. Rotor / number p. rotor		76								
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5094								
Radius / radius	mm	225								
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-11								

1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)

2) Angaben des Röhrchenherstellers beachten.

4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche


















5) 4255 nicht mit Deckel verschließbar



















1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)

2) Observe the tube manufacturer's instructions.

4) 4255-P: special surface treatment for highest hygienic requirements













5) 4255 cannot be closed with the lid
















4174	4579-A + 4255 / 4255-P 4)								
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\angle 90^\circ$  max. Laufzyklen /max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 								
	max. Laufzyklen /max. cycles: 30000 max. Beladung / max. load: 800 g (4500 RPM); 1000 g (4020 RPM), 1200 g (3670 RPM)								
	4435		4437		4438		4438 + 0726		
									
---	---	---	---	---	---	---	---		
									
Kapazität / capacity	ml	2,6 – 3,4	4,9	1,6 - 5	4 – 7	15	25	30	25
Maße / dimensions $\varnothing$ x L	mm	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	25 x 90	25 x 110	24 x 100
Anzahl p. Rotor / number p. rotor		84	84	84	84	48	28	28	28
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	<sup>2)</sup>	4935	4935	4935	4935	5207	5026	5026	4845
Radius / radius	mm	218	218	218	218	230	222	222	214
 g (97%)	sec	125							
 g	sec	197							
Temperatur / temperature	$^\circ\text{C}$ <sup>1)</sup>	-11							

4174	4579-A + 4255 / 4255-P 4)									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\angle 90^\circ$  max. Laufzyklen /max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 									
	max. Laufzyklen /max. cycles: 30000 max. Beladung / max. load: 800 g (4500 RPM); 1000 g (4020 RPM), 1200 g (3670 RPM)									
	4439		4440		4441		4442		4443	
										
---	Falcon®	Falcon®	---	---	---	---	5127			
										
Kapazität / capacity	ml	50	225	175	50	100	290	250		
Maße / dimensions $\varnothing$ x L	mm	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 137	61 x 122		
Anzahl p. Rotor / number p. rotor		16	4	4	20	8	4	4		
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500		
RZB / RCF	<sup>2)</sup>	4890	5501	5501	5207	4867	5320	5320		
Radius / radius	mm	216	243	243	230	215	235	235		
 g (97%)	sec	125								
 g	sec	197								
Temperatur / temperature	$^\circ\text{C}$ <sup>1)</sup>	-11								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 3) Gefäß nur belastbar bis RZB 700
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 5) 4255 nicht mit Deckel verschließbar
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 3) tube will not stand RCF values exceeding 700
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.













4174	4572				4547-B + 5621				
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    max. Laufzyklen / max. cycles: 30000 (4500RPM) 60000 (4000 RPM) 80000 (3500 RPM)									
	max. Laufzyklen / max. cycles: 19000 (4500 RPM), 45000 (3300 RPM) 90000 (2700 RPM)				max. Laufzyklen / max. cycles: 30000				
	max. Beladung / max. load: 1440 g				max. Beladung / max. load: 1200 g				
	---				---				
	4493		---		---		SK 32.08	SK 31.12-2	---
		---		---				---	---
---		---		---		Kartusche	Kartusche	---	---
				---				---	---
Kapazität / capacity	ml	1 – 5	4 - 7	---	---	---	---	---	---
Maße / dimensions	∅ x L	mm	13 x 75	13 x 100	---	---	43 x 193	43 x 114,5	---
Anzahl p. Rotor / number p. rotor		224		---	---	4	8	---	---
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	4718		---	---	5569	5524	---	---
Radius / radius	mm	208		---	---	246	244	---	---
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	3				-9			














4174	4547-B + 5621									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)										
	max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
			4449		4430		4432			
										
			0512	4239	Corning	Corning	Nalgene	Nunc	---	---
		 <sup>8)</sup>	 <sup>8)</sup>							
Kapazität / capacity	ml	750	1000	500	250	175	200	1,5	2,0	
Maße / dimensions	∅ x L	mm	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 139,2	60 x 130	11 x 38	11 x 38
Anzahl p. Rotor / number p. rotor		4	4	4	4	4	4	168	168	
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	
RZB / RCF	<sup>2)</sup>	5592	5592	5592	5592	5592	5592	5094	5094	
Radius / radius	mm	247	247	247	247	247	247	225	225	
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-9								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

































4174	4547-B + 5621									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\sphericalangle$ 90°  max. Laufzyklen /max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
	4433				4434					
										
	---	---	---	---	---	---	---	---		
	 2)	 2)			 2)		---	---		
Kapazität / capacity	ml	5	7	2,7 - 3	4,5 - 5	9	10	---	---	
Maße / dimensions	Ø x L	mm	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	16 x 80	---	---
Anzahl p. Rotor / number p. rotor		120				76			---	
Drehzahl / speed	RPM	4500							---	
RZB / RCF	<sup>2)</sup>	5026				5184			---	
Radius / radius	mm	222				229			---	
 9 (97%)	sec	125							---	
 9	sec	197							---	
Temperatur / temperature	°C <sup>1)</sup>	-9							---	

4174	4547-B + 5621									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\sphericalangle$ 90°  max. Laufzyklen /max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
	4434									
										
	---	---	---	---	---	---	---	---		
	 2)							 2)		
Kapazität / capacity	ml	8	4 - 5,5	7,5 - 8,2	9 - 10	4 - 7	8	5 - 10	12	
Maße / dimensions	Ø x L	mm	16 x 81	15 x 75	15 x 92	16 x 92	16 x 75	16 x 125	16 x 100	16 x 101
Anzahl p. Rotor / number p. rotor		76								
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5184								
Radius / radius	mm	229								
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-9								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.
















- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.
















4174	4547-B + 5621									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
	4434		4435				4437	4438 + 0726		
										
	---		---		---		---	---		
										
Kapazität / capacity	ml	10	15	2,6 – 2,9	4,9	1 - 5	4 - 7	15	25	
Maße / dimensions	Ø x L	mm	15 x 102	17 x 100	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	24 x 100
Anzahl p. Rotor / number p. rotor		76	76	84	84	84	84	48	28	
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5184	5184	5026	5026	5026	5026	5298	4913	
Radius / radius	mm	229	229	222	222	222	222	234	217	
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-9								

4174	4547-B + 5621									
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
	4438		4439		4440		4441	4442	4443	
										
	---		---		---		---	---	5127	
										
Kapazität / capacity	ml	25	30	50	225	175	50	100	250	
Maße / dimensions	Ø x L	mm	25 x 90	25 x 110	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 122
Anzahl p. Rotor / number p. rotor		28	28	16	4	4	20	8	4	
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5117	5117	4981	5592	5592	5298	4958	5411	
Radius / radius	mm	226	226	220	247	247	234	219	239	
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	-9								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4174	4547-B + 5621				4523-A						
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000				 max. Laufzyklen / max. cycles: 10000 (4500 RPM) 20000 (4000 RPM) 30000 (3500 RPM)						
	max. Beladung / max. load: 1200 g				max. Beladung / max. load: 1200 g						
	4443	---	---	SK 61.98	4548	4516-A		---			
		---	---					---			
	---	---	---	---	4-fach / 4-times	3-fach / 3-times	4-fach / 4-times	1-fach / 1-times	4589-A	1-fach / 1-times	
 8)	---	---		 13)	 13)	 12)	 13)				
Kapazität / capacity	ml	290	---	---	50	500	450	500	750	---	1000
Maße / dimensions	$\varnothing \times L$	mm	62 x 137	---	---	29 x 115	---	---	---	---	---
Anzahl p. Rotor / number p. rotor		4	---	---	20	4	---	4	---	---	4
Drehzahl / speed	RPM	4500									
RZB / RCF	<sup>2)</sup>	5411	---	---	5094	5524	---	5660	---	---	5705
Radius / radius	mm	239	---	---	225	244	---	250	---	---	252
 9 (97%)	sec	125									
 9	sec	197									
Temperatur / temperature	$^\circ\text{C}^{1)}$	-9							3		

4174	4523-A		---	SK 06.07 + 5629							
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 10000 (4500 RPM) 20000 (4000 RPM) 30000 (3500 RPM)		---	 max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)							
	max. Beladung / max. load: 1200 g		---	max. Beladung / max. load: 700 g							
	4508	SK 03.18	---	4626							
	 14)		---								
	4-fach / 4-times	1-fach / 1-times	4589-A	1-fach / 1-times	---	MTP	MTP	MS	CP		
 13)	 12)		---	---							
Kapazität / capacity	ml	500	750	---	300	---	---	---	---		
Maße / dimensions	TxBxH / DxVxH	mm	---	---	---	86x128x15	86x128x17,5	86x128x46	86x128x22		
Anzahl p. Rotor / number p. rotor		4		4	---	24	20	4	16		
Drehzahl / speed	RPM	4500									
RZB / RCF	<sup>2)</sup>	5660		5456	---	4324					
Radius / radius	mm	250		241	---	191					
 9 (97%)	sec	125									
 9	sec	197									
Temperatur / temperature	$^\circ\text{C}^{1)}$	3		---	---	---					

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 12) Ausgleichseinsatz. Handhabung siehe in Kapitel "Beladen des Rotors".
- 13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".
- 14) Einsatz mit Schlitz für Bänder zur Fixierung der Blutbeutel.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 12) Compensation insert. For handling see chapter "Loading the rotor".
- 13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".
- 14) Insert with slots for straps for fixing the blood bags












MTP Mikrotiterplatte /  
Microtitre plate
















CP Kulturplatte /  
Culture plate

DWP Deep Well Platte /  
Deep well plate

MS Micronic System /  
Micronic system

QP Filterplatte /  
Filter plate

4174		SK 06.07 + 5629						
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		 max. Laufzyklen / max. cycles: 30000 (4500 RPM), 60000 (4000 RPM), 80000 (3500 RPM)						
		max. Beladung / max. load: 700 g						
		4626		4626 + 1485		---	---	---
						---	---	---
DWP		Microtest-platten / plate Terasaki	QP	PCR-Platte 96-fach / PCR-plate 96 times	PCR-Strips	---	---	---
						---	---	---
Kapazität / capacity	ml	---	---	---	0,2	---	---	---
Maße / dimensions	mm	86x128x44,5	59x84x11	86x128x 83	82x124x20	---	---	---
Anzahl p. Rotor / number p. rotor		4	8	4	4	48	---	---
Drehzahl / speed	RPM	4500						
RZB / RCF	<sup>2)</sup>	4324						
Radius / radius	mm	191						
 9 (97%)	sec	125						
 9	sec	197						
Temperatur / temperature	°C <sup>1)</sup>	---						

4176		5052 + 5057								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		 max. Laufzyklen / max. cycles: 20000 (4500 RPM), 30000 (4000 RPM), 40000 (3500 RPM)								
		max. Beladung / max. load: 1100 g								
		mit Bioabdichtung / with bio-containment <sup>9)</sup>								
		4213		4213-93		4214				
										
---	---	---	---	---	---	---	---			
 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>	 <sup>2)</sup>			
Kapazität / capacity	ml	6	7	4,5 - 5	5	2,7 - 3	15	15	5 - 10	
Maße / dimensions	Ø x L	mm	12 x 82	12 x 100	11 x 92	12 x 75	11 x 66	17 x 100	17 x 100	16 x 100
Anzahl p. Rotor / number p. rotor		288			288		180			
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5818								
Radius / radius	mm	257								
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	14								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.

- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.

- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".















MTP Mikrotiterplatte /  
Microtitre plate













CP Kulturplatte /  
Culture plate

DWP Deep Well Platte /  
Deep well plate

MS Micronic System /  
Micronic system














QP Filterplatte /  
Filter plate




















4176		5052 + 5057									
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		 max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)									
		max. Beladung / max. load: 1100 g mit Bioabdichtung / with bio-containment <sup>9)</sup>									
		4214		4214-93		4215		4216		4218	
											
---		---		---		---		---			
											
Kapazität / capacity	ml	10	7,5 – 8,2	4 - 5,5	4 - 7	25	30	50	100		
Maße / dimensions $\varnothing$ x L	mm	15 x 102	15 x 92	15 x 75	16 x 75	24 x 100	25 x 110	34 x 100	44 x 100		
Anzahl p. Rotor / number p. rotor		180		180		66		36			
Drehzahl / speed	RPM	4500									
RZB / RCF	<sup>2)</sup>	5818									
Radius / radius	mm	257									
 9 (97%)	sec	125									
 9	sec	197									
Temperatur / temperature	°C <sup>1)</sup>	14									

4176		5052 + 5057									
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		 max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)									
		max. Beladung / max. load: 1100 g mit Bioabdichtung / with bio-containment <sup>9)</sup>									
		4220		4222		4222-93		---		---	
											
---		---		---		---		---			
											
Kapazität / capacity	ml	9 - 10	12	7	4 - 7	5	2,6 – 2,9	1 – 5	---		
Maße / dimensions $\varnothing$ x L	mm	16 x 92	17 x 100	12 x 100	13 x 100	12 x 75	13 x 65	13 x 75	---		
Anzahl p. Rotor / number p. rotor		96		180		180		---			
Drehzahl / speed	RPM	4500									
RZB / RCF	<sup>2)</sup>	5818									
Radius / radius	mm	257									
 9 (97%)	sec	125									
 9	sec	197									
Temperatur / temperature	°C <sup>1)</sup>	14									

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.














- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".















4176		5052 + 5057								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   $\sphericalangle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		 max. Laufzyklen / max. cycles: 20000 (4500 PM) 30000 (4000 RPM) 40000 (3500 RPM)								
		max. Beladung / max. load: 1100 g mit Bioabdichtung / with bio-containment <sup>9)</sup>								
		---	4223			---	4224		4225	
		---				---				
---	 <sub>2)</sub>		 <sub>2)</sub>	---	 <sub>2)</sub>			---		
Kapazität / capacity	ml	---	9	8	12	---	4	1,5	2	
Maße / dimensions $\varnothing \times L$	mm	---	14 x 100	16 x 125	16 x 101	---	10 x 88	11 x 38	11 x 38	
Anzahl p. Rotor / number p. rotor		---	150	150	150	150	378	504		
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5818			---		5818		5909	
Radius / radius	mm	---	257			---		257		261
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	14								

4176		5052 + 5057													
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   $\sphericalangle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		 max. Laufzyklen / max. cycles: 20000 (4500 PM) 30000 (4000 RPM) 40000 (3500 RPM)													
		max. Beladung / max. load: 1100 g mit Bioabdichtung / with bio-containment <sup>9)</sup>													
		4226	4232	4249	SK 09.00	4238		4241	4245-A						
															
Microtainer		---	---	---	5127		---	---	---						
								 <sub>8)</sub>		 <sub>8)</sub>					
Kapazität / capacity	ml	0,8	15	50	14	250	290	25	50						
Maße / dimensions $\varnothing \times L$	mm	8 x 45	17 x 120	29 x 115	16,5 x 106	62 x 122	62 x 137	25 x 90	29 x 115						
Anzahl p. Rotor / number p. rotor		468	138	36	138	6		48	48						
Drehzahl / speed	RPM	4500													
RZB / RCF	<sup>2)</sup>	5818	5999	5909	5954	5818		5818	5999						
Radius / radius	mm	257	265	261	263	257		257	265						
 9 (97%)	sec	125													
 9	sec	197													
Temperatur / temperature	°C <sup>1)</sup>	14													

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.


















- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".











4176		5052 + 5057				5052		---								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)										---						
		max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)						---		---						
		max. Beladung / max. load: 1100 g mit Bioabdichtung / with bio-containment <sup>9)</sup>						---		---						
		4258		4258 + 4449		6322		---		---		---				
								---		---		---				
0512		0554		4234-A		Corning		Corning		---						
										---						
Kapazität / capacity		ml		750		650		750		500		250				
Maße / dimensions		Ø x L		mm		97 x 152		96 x 139		96 x 135		96 x 147		60 x 162		
Anzahl p. Rotor / number p. rotor				6		6		6		6			---			
Drehzahl / speed		RPM		4500		4500		4500		4500			---			
RZB / RCF				5999		5999		5999		5818			---			
Radius / radius		mm		265		265		265		257			---			
		9 (97%)		sec		125		125		125			---			
		9		sec		197		197		197			---			
Temperatur / temperature		°C <sup>1)</sup>		14		14		14		14			---			

4176		4522-A																			
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)																					
		max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)																			
		max. Beladung / max. load: 1100 g ---																			
		4213		4213-93			4214			---											
										---											
---		---		---		---		---		---		---									
																					
Kapazität / capacity		ml		6		7		4,5 - 5		5		2,7 - 3		10		7,5 - 8,2		5 - 10			
Maße / dimensions		Ø x L		mm		12 x 82		12 x 100		11 x 92		12 x 75		11 x 66		15 x 102		15 x 92		16 x 100	
Anzahl p. Rotor / number p. rotor				288		288		288		288		180		180		180		180			
Drehzahl / speed		RPM		4500		4500		4500		4500		4500		4500		4500		4500			
RZB / RCF				5999		5999		5999		5999		5818		5818		5818		5818			
Radius / radius		mm		265		265		265		265		257		257		257		257			
		9 (97%)		sec		125		125		125		125		125		125		125			
		9		sec		197		197		197		197		197		197		197			
Temperatur / temperature		°C <sup>1)</sup>		14		14		14		14		14		14		14		14			

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 3) Gefäß nur belastbar bis RZB 700
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 3) tube will not stand RCF values exceeding 700
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".














4176		4522-A																	
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		  max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)  max. Beladung / max. load: 1100 g ---																	
		4214		4214-93		4215		4216		4218									
		 ---  2)  2)		 ---  		 ---  2) 		 ---  2)		 ---  2)									
Kapazität / capacity	ml	15		4 - 5,5		4 - 7		25		30		50		100					
Maße / dimensions	∅ x L	mm		17 x 100		17 x 100		15 x 75		16 x 75		24 x 100		25 x 110		34 x 100		44 x 100	
Anzahl p. Rotor / number p. rotor		180		180		180		66		36		24							
Drehzahl / speed	RPM	4500																	
RZB / RCF	2)	5818																	
Radius / radius	mm	257																	
 9 (97%)	sec	125																	
 9	sec	197																	
Temperatur / temperature	°C 1)	14																	




















4176		4522-A																	
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)		  max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)  max. Beladung / max. load: 1100 g ---																	
		4220		4222		4222-93		---		---									
		 --- 		 --- 		 --- 		---		---									
Kapazität / capacity	ml	9 - 10		12		7		4 - 7		5		2,6 - 2,9		1 - 5		---			
Maße / dimensions	∅ x L	mm		16 x 92		17 x 100		12 x 100		13 x 100		12 x 75		13 x 65		13 x 75		---	
Anzahl p. Rotor / number p. rotor		96		180		180		180		180		---		---		---			
Drehzahl / speed	RPM	4500																	
RZB / RCF	2)	5818																	
Radius / radius	mm	257																	
 9 (97%)	sec	125																	
 9	sec	197																	
Temperatur / temperature	°C 1)	14																	

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.



















4176	4522-A								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)								
	max. Beladung / max. load: 1100 g								
	---	4223			---	4224		4225	
	---				---				
---				---					
Kapazität / capacity	ml	---	9	8	12	---	4	1,5	2
Maße / dimensions $\varnothing$ x L	mm	---	14 x 100	16 x 125	16 x 101	---	10 x 88	11 x 38	11 x 38
Anzahl p. Rotor / number p. rotor		---	150			---	378	504	
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	---	5818			---	5818	5909	
Radius / radius	mm	---	257			---	257	261	
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	$^\circ\text{C}^{1)}$	14							

4176	4522-A								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)								
	max. Beladung / max. load: 1100 g								
	4226	4232	4249	SK 09.00	4238	4241	4245-A		
									
Microtainer	---	---	---	5127	---	---	---		
									
Kapazität / capacity	ml	0,8	15	50	14	250	290	25	50
Maße / dimensions $\varnothing$ x L	mm	8 x 45	17 x 120	29 x 115	16,5 x 106	62 x 122	62 x 137	25 x 90	29 x 115
Anzahl p. Rotor / number p. rotor		468	138	36	138	6		48	48
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	5818	5999	5909	5954	5818		5818	5999
Radius / radius	mm	257	265	261	263	257		257	265
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	$^\circ\text{C}^{1)}$	14							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4176	4522-A					4524-A				
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	  max. Laufzyklen / max. cycles: 20000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)  max. Beladung / max. load: 1100 g  ---					  max. Laufzyklen / max. cycles: 15000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)  max. Beladung / max. load: 2160 g  ---				
	4258		4258 + 4449		6322	4529-AO <sup>10)</sup>	4529-AM <sup>10)</sup>	4529-AU <sup>10)</sup>		
										
	0512	0554	4234-A	Corning	Corning	---	---	---		
	 8)		 8)		 8)		 13)			
	Kapazität / capacity	ml		750	650	750	500	250	500	
	Maße / dimensions	Ø x L		mm	97 x 152	96 x 139	96 x 135	96 x 147	60 x 162	---
	Anzahl p. Rotor / number p. rotor	6		6	6	6	12			
	Drehzahl / speed	RPM		4500						
	RZB / RCF	2)		5999	5999	5818	6498			
Radius / radius	mm		265	265	257	287				
 9 (97%)	sec		125							
 9	sec		197							
Temperatur / temperature	°C <sup>1)</sup>		14				16			

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 10) Darf nur in Gehänge 4524-A und beidseitig beladen verwendet werden. Bei Verwendung der Haltestifte darf die max. RZB 1000 nicht überschritten werden.


















	Position der Haltestifte	Abstand der Haltestifte vom Einsatzboden (innen)
4529-AO	oben	199,5 mm
4529-AM	Mitte	182,0 mm
4529-AU	unten	164,5 mm

- 12) Ausgleichseinsatz. Handhabung siehe in Kapitel "Beladen des Rotors".
- 13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 10) May only be used in bucket 4524-A and with both sides loaded. If the holding pins are used the max. RCF of 1000 may not be exceeded.

	Position of the holding pins	Distance of the holding pins from the bottom of the inserts (inside)
4529-AO	upper section	199.5 mm
4529-AM	middle section	182.0 mm
4529-AU	lower section	164.5 mm

- 12) Compensation insert. For handling see chapter "Loading the rotor".
- 13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".

4176	4524-A		4546-A		4591-A		4595-C						
<p><b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b></p>  <p>max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)</p>							 max. 3500 RPM zulässig						
	max. Laufzyklen / max. cycles: 15000 (4500-RPM) 30000 (4000 RPM) 40000 (3500 RPM)		max. Laufzyklen / max. cycle: 30000 (4500 RPM) 35000 (4000 RPM) 40000 (3500 RPM)		max. Laufzyklen / max. cycles: 15000 (4500 RPM) 30000 (4000 RPM) 40000 (3500 RPM)		max. Laufzyklen / max. cycles: 4500 (3500 RPM) 15000 (3000 RPM)						
	max. Beladung / max. load: 2160 g		max. Beladung / max. load: 2200 g		max. Beladung / max. load: 2160 g		max. Beladung / max. load: 2570 g						
	---		---		---		---						
	4592-B		4559-A		4598-A	4592-B		---	4596-A				
								---					
	4-fach / 4-times		4-fach / 4-times		3-fach / 3-times	4-fach / 4-times		0550	---				
	 13)		 12)		 13)	 12)		 8)	---				
	Kapazität / capacity	ml	500	750	---	450	---	450	500	750	---	2000	max. 1600
	Maße / dimensions	Ø x L	mm	---	---	---	---	---	---	---	---	180x150x100	---
Anzahl p. Rotor / number p. rotor		12		12		12	12		6	6			
Drehzahl / speed	RPM	4500		4500		4500	4500		3500	3500			
RZB / RCF		6498		6271		6498	6498		3848	3821			
Radius / radius	mm	287		277		287	287		281	279			
 9 (97%)	sec			125					95				
 9	sec			197					131				
Temperatur / temperature	°C <sup>1)</sup>			16								-3	

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 3) Gefäß nur belastbar bis RZB 700
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 10) Darf nur in Gehänge 4524-A und beidseitig beladen verwendet werden. Bei Verwendung der Haltestifte darf die max. RZB 1000 nicht überschritten werden.















	Position der Haltestifte	Abstand der Haltestifte vom Einsatzboden (innen)
4529-AO	oben	199,5 mm
4529-AM	Mitte	182,0 mm
4529-AU	unten	164,5 mm

















- 12) Ausgleichseinsatz. Handhabung siehe in Kapitel "Beladen des Rotors".
- 13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 3) tube will not stand RCF values exceeding 700
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 10) May only be used in bucket 4524-A and with both sides loaded. If the holding pins are used the max. RCF of 1000 may not be exceeded.

	Position of the holding pins	Distance of the holding pins from the bottom of the inserts (inside)
4529-AO	upper section	199.5 mm
4529-AM	middle section	182.0 mm
4529-AU	lower section	164.5 mm













- 12) Compensation insert. For handling see chapter "Loading the rotor".
- 13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".

















4176	4579-A	4579-A + 4255 / 4255-P <sup>4)</sup>								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)										
	max. Laufzyklen / max. cycles: 30000									
	---	max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RKM)								
	---	---	---	4449	4430			---		
	---	---	---					---		
	4255 / 4255-P <sup>4)</sup>	0512	4239	Corning	Corning	Nalgene®	Nunc®	---		
								---		
Kapazität / capacity	ml	1000	750	1000	500	250	175	200	---	
Maße / dimensions	∅ x L	mm	98 x 138	97 x 152	98 x 175	96 x 147	60 x 162	61,5 x 139,2	60 x 130	---
Anzahl p. Rotor / number p. rotor		6	6	6	6	6	6	6	---	
Drehzahl / speed	RPM	4500							---	
RZB / RCF	<sup>2)</sup>	6294							---	
Radius / radius	mm	278							---	
 9 (97%)	sec	125							---	
 9	sec	197							---	
Temperatur / temperature	°C <sup>1)</sup>	4							---	

4176	4579-A	4255 / 4255-P <sup>4)</sup>								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)										
	max. Laufzyklen / max. cycles: 30000									
	max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RPM)									
	4432	4433			4434					
										
---	---	---	---	---	---	---	---	---	---	
							---			
Kapazität / capacity	ml	1,5	2,0	5	7	2,7 - 3	4,5 - 5	9	---	15
Maße / dimensions	∅ x L	mm	11 x 38	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	---	17 x 100
Anzahl p. Rotor / number p. rotor		252		180			114			
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5796		5750			5886			
Radius / radius	mm	256		254			260			
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	4								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 5) 4255 nicht mit Deckel verschließbar
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 11) Maximale Beladung 800g. Bei einer Beladung über 800g muss die Drehzahl reduziert werden, siehe Beschriftung auf dem Becher. Berechnung der reduzierten Drehzahl siehe Kapitel "Zentrifugation von Stoffen oder Stoffgemischen mit einer höheren Dichte als 1,2 kg/dm<sup>3</sup>".




















- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 11) Maximum load 800g. With a load higher than 800g the speed has to be reduced, see label on the bucket. Calculation of the reduced speed see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>".








4176	4579-A + 4255 / 4255-P <sup>4)</sup>									
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 +									
	max. Laufzyklen / max. cycles: 30000									
	max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RPKM)									
	---									
	4434 									
	---	---	---	---	---	---	---	---		
	 5)		 2)					---		
Kapazität / capacity	ml	10	8	4 - 5,5	9 - 10	4 - 7	5 - 10	---		
Maße / dimensions	Ø x L	mm	15 x 102	16 x 80	16 x 81	15 x 75	16 x 92	16 x 75	16 x 100	---
Anzahl p. Rotor / number p. rotor		114							---	
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5886								
Radius / radius	mm	260								
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	4								

4176	4579-A + 4255 / 4255-P <sup>4)</sup>									
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   $\sphericalangle$ 90°  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 +									
	max. Laufzyklen / max. cycles: 30000									
	max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RPM)									
	---									
		4435		4437		4438		4438 + 0726		
										
	---	---	---	---	---	---	---	---		
								 2)		
Kapazität / capacity	ml	2,6 - 2,9	4,9	1 - 5	4 - 7	15	25	30	25	
Maße / dimensions	Ø x L	mm	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	25 x 90	25 x 110	24 x 100
Anzahl p. Rotor / number p. rotor		126			72		42		42	
Drehzahl / speed	RPM	4500								
RZB / RCF	<sup>2)</sup>	5750			6022		5818		5615	
Radius / radius	mm	254			266		257		248	
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	°C <sup>1)</sup>	4								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 5) 4255 nicht mit Deckel verschließbar
















- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid













4176	4579-A + 4255 / 4255-P <sup>4)</sup>							
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 + max. Laufzyklen / max. cycles: 30000							
	max. Beladung / max. load: 800g (4500 RPM), 1000g (4020 RPM), 1200g (3670 RPM)							
	---							
	4439	4440	4441	4442	4443	5127		SK 61.98
								
---	Falcon®	Falcon®	---	---	---	5127	---	
 2)	 5)	 5)	 5)	 2)	 5), 8)	 8)	 5)	
Kapazität / capacity ml	50	225	175	50	100	290	250	50
Maße / dimensions Ø x L mm	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 137	61 x 122	29 x 115
Anzahl p. Rotor / number p. rotor	24	6	30	12	6	6	6	30
Drehzahl / speed RPM	4500							
RZB / RCF <sup>2)</sup>	5705	6294	5999	5683	6113	6113	6113	5796
Radius / radius mm	252	278	265	251	270	270	270	256
 9 (97%) sec	125							
 9 sec	197							
Temperatur / temperature °C <sup>1)</sup>	4							

4176	4572	---						
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 --- max. Laufzyklen / max. cycles: 19000 (4500 RPM), 45000 (3300 RPM) 90000 (2700 RPM) --- max. Beladung / max. load: 1440 g --- ---							
	---							
	---							
	4493							
 --- --- --- ---								
---								
 1-5	 4-7	---	---	---	---	---	---	
Kapazität / capacity ml	1 – 5	4 – 7	---	---	---	---	---	
Maße / dimensions Ø x L mm	13 x 75	13 x 100	---	---	---	---	---	
Anzahl p. Rotor / number p. rotor	336							
Drehzahl / speed RPM	4500							
RZB / RCF <sup>2)</sup>	5522							
Radius / radius mm	244							
 9 (97%) sec	125							
 9 sec	197							
Temperatur / temperature °C <sup>1)</sup>	10							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur (nur bei Kühlzentrifuge)
- 2) Angaben des Röhrchenherstellers beachten.
- 3) Gefäß nur belastbar bis RZB 700
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 5) 4255 nicht mit Deckel verschließbar
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.







- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Observe the tube manufacturer's instructions.
- 3) tube will not stand RCF values exceeding 700
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.













4176	4547-B + 5621									
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   $\sphericalangle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
	---		4449		4430		4432			
	---									
	0512	4239	Corning	Corning	Nalgene	Nunc	---	---		
	 8)	 8)								
	Kapazität / capacity	ml	750	1000	500	250	175	200	1,5	2,0
	Maße / dimensions $\varnothing \times L$	mm	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 139,2	60 x 130	11 x 38	11 x 38
	Anzahl p. Rotor / number p. rotor		6	6	6	6	6	6	6	252
	Drehzahl / speed	RPM	4500							
	RZB / RCF	<sup>2)</sup>	6384	6384	6384	6384		5886		
Radius / radius	mm	282	282	282	282		260			
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	$^\circ\text{C}^{1)}$	9								

4176	4547-B + 5621									
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   $\sphericalangle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g									
	4433				4434					
										
	---	---	---	---	---	---	---	---		
	 2)	 2)			 2)	---		---		
	Kapazität / capacity	ml	5	7	2,7 - 3	4,5 - 5	9	---	10	---
	Maße / dimensions $\varnothing \times L$	mm	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	---	16 x 80	---
	Anzahl p. Rotor / number p. rotor		180				114		---	
	Drehzahl / speed	RPM	4500							
	RZB / RCF	<sup>2)</sup>	5841				5977		---	
Radius / radius	mm	258				264		---		
 9 (97%)	sec	125								
 9	sec	197								
Temperatur / temperature	$^\circ\text{C}^{1)}$	9								

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 4) 4255-P: spezielle Oberflächenbehandlung für höchste hygienische Ansprüche
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.












4176	4547-B + 5621								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 - 5RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g								
	4434								
									
									
Kapazität / capacity	ml	8	4 - 5,5	7,5 - 8,2	9 - 10	4 - 7	8	5 - 10	12
Maße / dimensions $\varnothing \times L$	mm	16 x 81	15 x 75	15 x 92	16 x 92	16 x 75	16 x 125	16 x 100	16 x 101
Anzahl p. Rotor / number p. rotor		114							
Drehzahl / speed	RPM	4500							
RZB / RCF	<sup>2)</sup>	5977							
Radius / radius	mm	264							
 9 (97%)	sec	125							
 9	sec	197							
Temperatur / temperature	°C <sup>1)</sup>	9							



















4176	4547-B + 5621										
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g										
	4434		4435				4437		4438 + 0726		
											
											
Kapazität / capacity	ml	10	15	2,6 - 2,9	4,9	1 - 5	4 - 7	15	25		
Maße / dimensions $\varnothing \times L$	mm	15 x 102	17 x 100	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	24 x 100		
Anzahl p. Rotor / number p. rotor		114				126					
Drehzahl / speed	RPM	4500									
RZB / RCF	<sup>2)</sup>	5977				5841		6090		5728	
Radius / radius	mm	264				258		269		253	
 9 (97%)	sec	125									
 9	sec	197									
Temperatur / temperature	°C <sup>1)</sup>	9									

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.












4176	4547-B + 5621											
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g											
	4438		4439		4440		4441		4442		4443	
												
	---		---		---		---		---		5127	
Kapazität / capacity	ml	25	30	50	225	175	50	100	250			
Maße / dimensions	Ø x L	mm	25 x 90	25 x 110	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 122		
Anzahl p. Rotor / number p. rotor		42		24		6		30		6		
Drehzahl / speed	RPM	4500										
RZB / RCF	<sup>2)</sup>	5909		5773		6384		6090		6203		
Radius / radius	mm	261		255		282		269		274		
 9 (97%)	sec	125										
 9	sec	197										
Temperatur / temperature	°C <sup>1)</sup>	9										

4176	4547-B + 5621						4523-A										
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)	 max. Laufzyklen / max. cycles: 30000 max. Beladung / max. load: 1200 g						 max. Laufzyklen / max. cycles: 10000 (4500 RPM) 20000 (4000 RPM) 30000 (3500 RPM) max. Beladung / max. load: 1200 g										
	4443		---		SK 20.16		SK 61.98		4548		4516-A		---				
			---										---				
	---		---		Flasche / bottle		---		4-fach / 4-times		3-fach / 3-times		4-fach / 4-times		1-fach / 1-times		
		---															
Kapazität / capacity	ml	290		---		375		50		500		450		---		1000	
Maße / dimensions	Ø x L	mm		62 x 137		---		65 x 145,5		29 x 115		---		---		---	
Anzahl p. Rotor / number p. rotor		6		---		---		30		6		6		6		6	
Drehzahl / speed	RPM	4500															
RZB / RCF	<sup>2)</sup>	6203		---		6294		5886		6316		6475		6520			
Radius / radius	mm	274		---		278		260		279		286		288			
 9 (97%)	sec	125															
 9	sec	197															
Temperatur / temperature	°C <sup>1)</sup>	9															

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 12) Ausgleichseinsatz. Handhabung siehe in Kapitel "Beladen des Rotors".
- 13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 12) Compensation insert. For handling see chapter "Loading the rotor".
- 13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".

4176	4523-A		---	---	---	---	---
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>    $\angle 90^\circ$  max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)			---	---	---	---	---
	max. Laufzyklen / max. cycles: 10000 (4500 RPM) 20000 (4000 RPM) 30000 (3500 RPM)		---	---	---	---	---
	max. Beladung / max. load: 1200 g ---		---	---	---	---	---
	4508      SK 03.18		---	---	---	---	---
	 14)			---	---	---	---
4-fach / 4-times    1-fach / 1-time    4589-A    1-fach		---	---	---	---	---	
 13)  12)  13)		---	---	---	---	---	
Kapazität / capacity      ml	500	750	---	300	---	---	---
Maße / dimensions $\varnothing \times L$ mm	---	---	---	---	---	---	---
Anzahl p. Rotor / number p. rotor	6		---	6	---	---	---
Drehzahl / speed                  RPM	4500		---	---	---	---	
RZB / RCF                              2)	6475		---	6271	---	---	
Radius / radius                        mm	286		---	277	---	---	
 9 (97%)                              sec	125		---	---	---	---	
 9    sec	197		---	---	---	---	
Temperatur / temperature        °C 1)	9		---	---	---	---	

1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur

2) Angaben des Röhrchenherstellers beachten.

12) Ausgleichseinsatz. Handhabung siehe in Kapitel "Beladen des Rotors".

13) Ausgleichsgewichte Set (4566) erhältlich. Handhabung siehe in Kapitel "Beladen des Rotors".

14) Einsatz mit Schlitz für Bänder zur Fixierung der Blutbeutel.













1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature

2) Observe the tube manufacturer's instructions.

12) Compensation insert. For handling see chapter "Loading the rotor".

13) Balance weights set (4566) available. For handling see chapter "Loading the rotor".

14) Insert with slots for straps for fixing the blood bags

4176		SK 06.07 + 5629								
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 RPM) 60000 (4000 RPM) 80000 (3500 RPM)										
		max. Laufzyklen / max. cycles: 30000 (4500 RPM), 60000 (4000 RPM), 80000 (3500 RPM)								
		max. Beladung / max. load: 700 g ---								
		<b>4626</b> 								
		MTP	MTP	MS	CP	DWP	Microtest- platten / plate Terasaki	QP	---	
									---	
Maße / dimensions TxBxH / DxWxH	mm	86x128x15	86x128x17,5	86x128x46	86x128x22	86x128x44,5	59x84x11	86x128x83	---	
Anzahl p. Rotor / number p. rotor		36	30	6	24	6	12	6	---	
Drehzahl / speed	RPM								4500	---
RZB / RCF	<sup>2)</sup>								5139	---
Radius / radius	mm								227	---
 9 (97%)	sec								125	---
 9	sec								197	---
Temperatur / temperature	°C <sup>1)</sup>								-3	---

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.








MTP Mikrotiterplatte /  
Microtitre plate

















CP Kulturplatte /  
Culture plate

DWP Deep Well Platte /  
Deep well plate

MS Micronic System /  
Micronic system




QP Filterplatte /  
Filter plate

















4176		SK 06.07 + 5629							
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>   max. Laufzyklen / max. cycles: 30000 (4500 - 4001 RPM) 60000 (4000 - 3501 RPM) 80000 (3500 - 50 RPM)									
		max. Laufzyklen / max. cycles: 30000 (4500 - 4001 RPM), 60000 (4000 - 3501 RPM), 80000 (3500 - 50 RPM)							
		max. Beladung / max. load: 700 g ---							
		<b>4626 + 1485</b>  							
		PCR-Platte 96-fach / PCR-plate 96 times	PCR-Strips	---	---	---	---	---	---
				---	---	---	---	---	---
Kapazität / capacity	ml	---	8 x 0,2	---	---	---	---	---	---
Maße / dimensions TxBxH / DxWxH	mm	82x124x20	---	---	---	---	---	---	---
Anzahl p. Rotor / number p. rotor		6	72	---	---	---	---	---	---
Drehzahl / speed	RPM	4500	---	---	---	---	---	---	---
RZB / RCF	<sup>2)</sup>	5139	---	---	---	---	---	---	---
Radius / radius	mm	227	---	---	---	---	---	---	---
 9 (97%)	sec	125	---	---	---	---	---	---	---
 9	sec	197	---	---	---	---	---	---	---
Temperatur / temperature	°C <sup>1)</sup>	- 3	---	---	---	---	---	---	---

4570		---							
<b>Winkelrotor 6-fach / Angle rotor 6-times</b>   max. Laufzyklen / max. cycles: 15000  Hülsen / reduction (6x inclusive) max. Laufzyklen: 15000 einsetzbar bis / usable until: 5 Jahre / year  mit Bioabdichtung / with bio-containment <sup>9)</sup>		---							
		---							
		---							
		---							
		---	5641	5642	5643	5644	---	---	---
		---					---	---	---
		5127	---	---	---	---	---	---/---	---
		 8)			 2)			 + 	 8)
Kapazität / capacity	ml	250	10	30	25	50	85	94	94
Maße / dimensions Ø x L	mm	61,5 x 122	16 x 80	26 x 95	24 x 100	29 x 107	38 x 106	38 x 110	38 x 102
Anzahl p. Rotor / number p. rotor		6	48	18	6	6	6	6	6
Drehzahl / speed	RPM	6000							
RZB / RCF	<sup>2)</sup>	5594	5353	5152	4830	4910	4910	4910	4910
Radius / radius	mm	139	133	128	120	122	122	122	122
 9 (97%)	sec	64							
 9	sec	69							
Temperatur / temperature	°C <sup>1)</sup>	- 16							

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrenherstellers beachten.
- 8) Bei Temperaturen über 40 °C und/oder geringer Befüllung der Gefäße können sich diese verformen.
- 9) Nach DIN EN 61010, Teil 2 – 020. Die Hinweise für Bio-Sicherheitssysteme in den Kapiteln "Sicherheitshinweise" und "Pflege und Wartung" beachten.







- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".







4570		---						
<b>Winkelrotor 6-fach / Angle rotor 6-times</b>  $\angle 25^\circ$	---							
	---							
max. Laufzyklen / max. cycles: 15000	---							
Hülsen / reduction (6x inclusive) max. Laufzyklen: 15000 einsetzbar bis / usable until: 5 Jahre / year mit Bioabdichtung / with bio-containment <sup>9)</sup>	---							
Kapazität / capacity ml	15	15	15	50	---	---	---	---
Maße / dimensions $\varnothing \times L$ mm	17 x 100	17 x 100	17 x 120	29 x 115	---	---	---	---
Anzahl p. Rotor / number p. rotor	42		30	6	---	---	---	---
Drehzahl / speed RPM	6000			---				
RZB / RCF <sup>2)</sup>	5313		5152	4830	---	---	---	---
Radius / radius mm	132		128	120	---	---	---	---
 9 (97%) sec	64			---				
 9 sec	69			---				
Temperatur / temperature $^\circ\text{C}$ <sup>1)</sup>	- 16			---				

SK 48.06-35		SK 48.06-36						---	
<b>Ausschwingrotor 4-fach / Swing out rotor 4-times</b>  $\angle 90^\circ$							---		
	max. Laufzyklen: 30000 2800 RPM max. load 1300g 1000 RPM max. load 2300g						---		
---		---						---	
	<b>SK 48.06-3</b>	<b>SK 48.06-4</b>	<b>SK 48.06-7</b>	<b>SK 48.06-11</b>	<b>SK 48.06-12</b>	<b>SK 48.06-13</b>	---	---	
							---	---	
	Kartusche	Kartusche	Kartusche	Kartusche	Kartusche	Kartusche	---	---	
max. Laufzyklen / max. cycles: 30000							---	---	
Kapazität / capacity ml	---	---	---	---	---	---	---	---	
Maße / dimensions $\varnothing \times L$ mm	18,5 x 91	25,5 x 118	42,5 x 186	18,65 x 95,6	25,6 x 122,4	26 x 180,4	---	---	
Anzahl p. Rotor / number p. rotor	28	16	8	28	12	12	---	---	
Drehzahl / speed RPM	2800						---		
RZB / RCF <sup>2)</sup>	2375	2375	2323	2375	2375	2375	---	---	
Radius / radius mm	271	271	265	271	271	271	---	---	
 9 (97%) sec	---						---		
 9 sec	---						---		
Temperatur / temperature $^\circ\text{C}$ <sup>1)</sup>	7						---		

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur
- 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Observe the tube manufacturer's instructions.

SK 08.09	SK 08.09-5	---		---		---	
<b>Ausschwingrotor 2-fach / Swing out rotor 2-times</b>  max. Laufzyklen / max. cycles: 100000		---	---	---	---	---	---
	max. Laufzyklen / max. cycles: 20000 max. Beladung / max. load: 1915g	---	---	---	---	---	---
	---	---	---	---	---	---	---
	<b>SK 08.09-4</b>	---	---	---	---	---	---
		---	---	---	---	---	---
	<b>STIWA-Rack</b>	---	---	---	---	---	---
		---	---	---	---	---	---
Kapazität / capacity	ml	---	---	---	---	---	---
Maße / dimensions $\varnothing$ x L	mm	22 x 60 x 220	---	---	---	---	---
Anzahl p. Rotor / number p. rotor		10	---	---	---	---	---
Drehzahl / speed	RPM	3520	---				
RZB / RCF	<sup>2)</sup>	2812	---	---	---	---	---
Radius / radius	mm	203	---	---	---	---	---
 9 (97%)	sec	82	---				
 9	sec	137	---				
Temperatur / temperature	$^{\circ}\text{C}$ <sup>1)</sup>	---	---				

SK 48.08-1	SK 48.08	---		---		---	
<b>Ausschwingrotor 6-fach / Swing out rotor 6-times</b>  max. Laufzyklen / max. cycles: 30000		---	---	---	---	---	---
	max. Laufzyklen / max. cycles: 50000 max. Beladung / max. load: 1130g	---	---	---	---	---	---
	---	---	---	---	---	---	---
	<b>SK 48.08-2</b>	---	---	---	---	---	---
		---	---	---	---	---	---
	<b>Nashua Kits</b>	---	---	---	---	---	---
		---	---	---	---	---	---
Kapazität / capacity	ml	---	---	---	---	---	---
Maße / dimensions $\varnothing$ x L	mm	86 x 128 x 44,5	---	---	---	---	---
Anzahl p. Rotor / number p. rotor		6	---	---	---	---	---
Drehzahl / speed	RPM	3400	---				
RZB / RCF	<sup>2)</sup>	3270	---	---	---	---	---
Radius / radius	mm	253	---	---	---	---	---
 9 (97%)	sec	---	---				
 9	sec	---	---				
Temperatur / temperature	$^{\circ}\text{C}$ <sup>1)</sup>	---	---				

- 1) Tiefste erreichbare Temperatur bei maximaler Drehzahl, 1 h Laufzeit und 20°C Raumtemperatur  
 2) Angaben des Röhrchenherstellers beachten.

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature  
 2) Observe the tube manufacturer's instructions.